


Tooth Extracting Forceps|(eng)



Incisors and Cuspids 1

XU-51-001

Tooth Ext Forceps
Fig.1,



Incisors and Bicusps 2

XU-51-002


Tooth Ext Forceps
Fig.2,



Lower Incisors and Canines 4

XU-51-004

Tooth Ext Forceps
Fig.4,



Bicusps 7

XU-51-007


Tooth Ext Forceps
Fig.7,



Lower Bicusps, Universal 13

XU-51-013

Tooth Ext Forceps
Fig.13,




Molars, Right 17

XU-51-017

Tooth Ext Forceps
Fig.17,

Tooth Extracting Forceps|(eng)



Molars, Left 18
XU-51-018
Tooth Ext Forceps
Fig.18,



Molars 18A
-51-018-01
Tooth Ext Forceps
Fig.18A,



Lower Molars, Universal 22
XU-51-022
Tooth Ext Forceps
Fig.22,



Lower Molars, Right 23
XU-51-023
Tooth Ext Forceps
Fig.23,



Lower Molars, Left 24
XU-51-024
Tooth Ext Forceps
Fig.24,



Upper roots and Incisors 29
XU-51-029
Tooth Ext Forceps
Fig.29,

Tooth Extracting Forceps|(eng)



Lower Roots 33
XU-51-033
Tooth Ext Forceps
Fig.33,



Lower Roots 01
XU-51-033-01
Tooth Ext Forceps
Fig.33 A,



Lower Roots 03
XU-51-033-03
Tooth Ext Forceps
Fig.33 C,



Lower Roots 06
XU-51-033-06
Tooth Ext Forceps
Fig.33 M,



Incisors and Bicuspid 34N
XU-51-034
Tooth Ext Forceps
Fig.34,



Incisors and Bicuspid 35N
XU-51-035
Tooth Ext Forceps
Fig.35,

Tooth Extracting Forceps|(eng)



XU-51-036
Tooth Ext Forceps
Fig.36,



Lower roots 46
XU-51-046
Tooth Ext Forceps
Fig.46,



Lower Roots 46
XU-51-046-01
Tooth Ext Forceps
Fig.46 L,



Lower Roots 46
XU-51-046-02
Tooth Ext Forceps
Fig.46 LX,



Very Fine Roots 49
XU-51-049
Tooth Ext Forceps
Fig.49,



Roots 51
-51-051
Tooth Ext Forceps
Fig.51,

Tooth Extracting Forceps|(eng)



Roots 51A
-51-051-01
Tooth Ext Forceps
Fig.51 A,



Roots 51 C
-51-051-03
Tooth Ext Forceps
Fig.51 C,



XU-51-057
Tooth Ext Forceps
Fig.57,



XU-51-057-1
Tooth Ext Forceps
Fig.57 A,



XU-51-057-02
Tooth Ext Forceps
Fig.57 B,



Third Molars 67
XU-51-067
Tooth Ext Forceps
Fig.67,

Tooth Extracting Forceps|(eng)



Third Molars 67 A
XU-51-067-01
Tooth Ext Forceps
Fig.67 A,



Third Molars Deep Gripping 67 N
XU-51-067-02
Tooth Ext Forceps
Fig.67 N,



Lower Molars 73
XU-51-073
Tooth Ext Forceps
Fig.73,



Lower Roots 74
XU-51-074
Tooth Ext Forceps
Fig.74,



Incisors and Roots 74N
XU-51-074-06
Tooth Ext Forceps
Fig.74 N,



Incisors and Roots ,Slender Jaw74XN
XU-51-074-07
Tooth Ext Forceps
Fig.74 XN,

Tooth Extracting Forceps|(eng)



Incisors and Bicuspids
XU-51-075
Tooth Ext Forceps
Fig.75,



Upper Roots 76
XU-51-076
Tooth Ext Forceps
Fig.76,



Upper Roots 76N
XU-51-076-07
Tooth Ext Forceps
Fig.76 N,



Third Molars 79
XU-51-079
Tooth Ext Forceps
Fig.79,



Third Molars 79A
XU-51-079-01
Tooth Ext Forceps
Fig.79 A,



Third Molars ,Deep Gripping 79N
XU-51-079-02
Tooth Ext Forceps
Fig.79 N,

Tooth Extracting Forceps|(eng)



Molars with carious or broken caps 86A

XU-51-086-01

Tooth Ext Forceps
Fig.86 A,



Molars with carious or broken caps 86B

XU-51-086-02

Tooth Ext Forceps
Fig.86 B,



Molars with carious or broken caps 86C

XU-51-086-03

Tooth Ext Forceps
Fig.86 C,



XU-51-097

Tooth Ext Forceps
Fig.97,



Very Fine Roots 133

XU-51-133

Tooth Ext Forceps
Fig.133,



Roots Fragments, Deep Gripping 146

XU-51-146

Tooth Ext Forceps
Fig.146,

Tooth Extracting Forceps|(eng)



Tooth Ext Forceps
XU-51-151
Tooth Ext Forceps
Fig.151,



Tooth Ext Forceps
XU-51-153
Fig.153,



Roots, with Parallel beak 233
XU-51-233
Tooth Ext Forceps
Fig.233,



Incisors and Bicuspid 1
XU-51-301
Tooth Ext Forceps
Fig.1,



Incisors and Bicuspid 2
XU-51-302
Tooth Ext Forceps
Fig.2,



Bicuspid 7
XU-51-307
Tooth Ext Forceps
Fig.7,

Tooth Extracting Forceps|(eng)



Bicuspid 13
XU-51-313
Tooth Ext Forceps
Fig.13,



Molars 17
-51-317
Tooth Ext Forceps
Fig.17,



Molars 18
-51-318
Tooth Ext Forceps
Fig.18,



Molars 22
-51-322
Tooth Ext Forceps
Fig.22,



Roots 33
-51-330
Tooth Ext Forceps
Fig.30,



Roots 33
-51-333
Tooth Ext Forceps
Fig.33,

Tooth Extracting Forceps|(eng)



Roots 51
-51-351
Tooth Ext Forceps
Fig.51,



Roots, Very Fine Beak 51A
XU-51-351-01
Tooth Ext Forceps
Fig.51 A,



Third Molars 67 A
XU-51-367-01
Tooth Ext Forceps
Fig.67 A,



Molars 73
-51-373
Tooth Ext Forceps
Fig.73,




Incisors and Roots 74
XU-51-374
Tooth Ext Forceps
Fig.74,



Incisors and roots , with fine beak 74N
XU-51-374-01
Tooth Ext Forceps
Fig.74 N,

Tooth Extracting Forceps|(eng)



Third Molars 79
XU-51-379
Tooth Ext Forceps
Fig.79,



Incisors and cuspids 1
XU-51-401
Tooth Ext Forceps
Fig.1 Diam,



Bicuspid 7
XU-51-407
Tooth Ext Forceps
Fig.7 Diam,



Incisors and cuspids 13
XU-51-413
Tooth Ext Forceps
Fig.13 Diam,



Molars , Right 17
XU-51-417
Tooth Ext Forceps
Fig.17 Diam,



Molars , left 18
XU-51-418
Tooth Ext Forceps
Fig.18 Diam,

Tooth Extracting Forceps|(eng)



Molars 22

-51-422

Tooth Ext Forceps
Fig.22 Diam,



Roots 33

-51-433

Tooth Ext Forceps
Fig.33 Diam,



Very Fine Roots 46

XU-51-446

Tooth Ext Forceps
Fig.46 Diam,



Very Fine Roots 49

XU-51-449

Tooth Ext Forceps
Fig.49 Diam,



Roots 51

-51-451

Tooth Ext Forceps
Fig.51 Diam,



Roots , With Fine Beak 51A

XU-51-451-01

Tooth Ext Forceps
Fig.51 A Diam,

Tooth Extracting Forceps|(eng)



Martin Ergo 67A
XU-51-467-01
Tooth Ext Forceps
Fig.67 A Diam,



Martin Ergo 74N
XU-51-474
Tooth Ext Forceps
Fig.74 N Diam,



Martin Ergo 79
XU-51-479
Tooth Ext Forceps
Fig.79 Diam,



Bicuspid 13S
XU-51-613-09
Tooth Ext Forceps
Fig.13 S,



Molars 22s
-51-622-09
Tooth Ext Forceps
Fig.22 S,

Martin Ergo ,Very Fine Roots 97
XU-51-497
Tooth Ext Forceps
Fig.97 Diam,

Tooth Extracting Forceps|(eng)



Incisors and Roots 30S

XU-51-630-09

Tooth Ext Forceps
Fig.30 S,



Roots 33S

-51-633-09

Tooth Ext Forceps
Fig.33 S,



Incisors and bicuspid 37

XU-51-637

Tooth Ext Forceps
Fig.37,



Incisors and Cuspid 38

XU-51-638

Tooth Ext Forceps
Fig.38,



Molars 39

XU-51-

Tooth Ext Forceps
Fig.39,



Molars 39A

XU-51-639-

Tooth Ext Forceps
Fig.39 A,

Tooth Extracting Forceps|(eng)



Molars , Left 39L
XU-51-639-05
Tooth Ext Forceps
Fig.39 L,



Molars , Right 39R
XU-51-639-08
Tooth Ext Forceps
Fig.39 R,



Molars 40
-51-640
Tooth Ext Forceps
Fig.40,



Roots 51S
-51-651-09
Tooth Ext Forceps
Fig.51 S,



Roots 76S
-51-676-09
Tooth Ext Forceps
Fig.76 S,



Molars 3
-51-703
Tooth Ext Forceps
Fig.3,

Tooth Extracting Forceps|(eng)



Incisors 5
-51-705
Tooth Ext Forceps
Fig.5,



Molars 6
-51-706
Tooth Ext Forceps
Fig.6,



Roots 7
-51-707
Tooth Ext Forceps
Fig.7,



Incisors and cuspids 137
XU-51-737
Tooth Ext Forceps
Fig.137,



Bicuspid 139
XU-51-739
Tooth Ext Forceps
Fig.139,



Klein , Roots 51S
XU-51-751-09
Tooth Ext Forceps
Fig.51 S,

Tooth Extracting Forceps|(eng)



Molars 17SK
XU-51-817-09
Tooth Ext Forceps
Fig.17 SK,



Molars 23SK
XU-51-823-09
Tooth Ext Forceps
Fig.23 SK,



Teeth and Roots ,Universal 150SK
XU-51-850-09
Tooth Ext Forceps
Fig.150 SK,



Teeth and Roots, Universal 151Sk
XU-51-851-09
Tooth Ext Forceps
Fig.151 SK,

Tooth Extracting Forceps|(amr)



Incisors and Bicuspid 1

VI-52-001

Tooth Ext Forceps Amr
Fig.1,



Henahan Incisors 1A

VI-52-001-01

Tooth Ext Forceps Amr
Fig.1 A,



Molars and third molars 10S

XU-52-010-09

Tooth Ext Forceps Amr
Fig.10 S,



Molars 15

-52-015

Tooth Ext Forceps Amr
Fig.15,



Molars 16

-52-016

Tooth Ext Forceps Amr
Fig.16,



Molars , Slender Jaw

XU-52-016-09

Tooth Ext Forceps Amr
Fig.16 S,

Tooth Extracting Forceps(amr)



Molars 17
-52-017
Tooth Ext Forceps Amr
Fig.17,



Harris Molars , left 18L
XU-52-018-05
Tooth Ext Forceps Amr
Fig.18 L,



Harris Molars , Right 18R
XU-52-018-08
Tooth Ext Forceps Amr
Fig.18 R,



Molars
-52-023
Tooth Ext Forceps Amr
Fig.23,



Parnly Molars and bicuspid 32
XU-52-032
Tooth Ext Forceps Amr
Fig.32,



Parnly Molars and Bicuspid , Roots 32A
XU-52-032-01
Tooth Ext Forceps Amr
Fig. 32 A,

Tooth Extracting Forceps|(amr)



Molars , Left 53L
XU-52-053-05
Tooth Ext Forceps Amr
Fig.53 L,



Molars , Right 53R
XU-52-053-08
Tooth Ext Forceps Amr
Fig.53 R,



**Teeth and Roots, Universal also for
Children 62**
XU-52-062
Tooth Ext Forceps Amr
Fig.62,



Incisors and Roots 65
XU-52-065
Tooth Ext Forceps Amr
Fig.65,



**Tomes Roots Fragments and Small Root
69**
XU-52-069
Tooth Ext Forceps Amr
Fig.69,



Nevius Molars , Left 88L
XU-52-088-05
Tooth Ext Forceps Amr
Fig.88 L,

Tooth Extracting Forceps[(amr)



Nevius Molars , Right 88R
XU-52-088-08

Tooth Ext Forceps Amr
Fig.88 R,



Cryer Bicuspid, Incisors and Roots 150
XU-52-150

Tooth Ext Forceps Amr
Fig.150,



Cryer Bicuspid, Incisors and Roots 150A
XU-52-150-01

Tooth Ext Forceps Amr
Fig.150 A,



Bicuspid and Incisors , Roots 150X
XU-52-150-02

Tooth Ext Forceps Amr
Fig.150 X,



Bicuspid and Incisors, Roots 150AS
XU-52-150-08

Tooth Ext Forceps Amr
Fig.150 AS,



Cryer Bicuspid and Incisors , Roots 150S
XU-52-150-09

Tooth Ext Forceps Amr
Fig.150 S,

Tooth Extracting Forceps|(amr)



Cryer Bicuspid, Incisors and Roots 151

XU-52-151

Tooth Ext Forceps Amr
Fig.151,



Cryer Bicuspid, Incisors and Roots 151A

XU-52-151-01

Tooth Ext Forceps Amr
Fig.151 A,



Cryer Bicuspid, Incisors and Roots 151X

XU-52-151-02

Tooth Ext Forceps Amr
Fig.151 X,



Bicuspid and Incisors, Roots 151AS

XU-52-151-08

Tooth Ext Forceps Amr
Fig.151 AS,



Cryer Teeth and Roots, Universal 151S

XU-52-151-09

Tooth Ext Forceps Amr
Fig.151 S,



Third Molars 222

XU-52-222

Tooth Ext Forceps Amr
Fig.222,

Tooth Extracting Forceps|(amr)



Root Fragments and Universal 501

XU-52-501

Tooth Ext Forceps Amr
Fig.501,



Witzel Root Fragments and Universal 502

XU-52-502

Tooth Ext Forceps Amr
Fig.502,

Root Elevators




Seldin 1L
XU-53-001-05
Root Elevators
Fig.1L,




Seldin 1R
XU-53-001-08
Root Elevator
Fig.1R,



Seldin 2
-53-002-00
Root Elevators
Fig.2,



Seldin
-53-003-00
Root Elevators
Fig.3,



Cryer
-53-003-01
Root Elevators
Fig.3 A,



Seldin 4L
XU-53-004-05
Root Elevators
Fig. 4 L,

Root Elevators



Seldin 4R
XU-53-004-08
Root Elevators
Fig. 4 R,



Crane Pick 6
XU-53-006-00
Root Elevators
Fig. 6,



Miller 1
-53-008-00
Root Elevators
Fig. 8,



Miller 2
-53-009-00
Root Elevators
Fig.9,



Miller-Apexo
XU-53-010-00
Root Elevators
Fig.10,



Miller-Apexo
XU-53-011-00
Root Elevators
Fig.11,

Root Elevators



Cryer
-53-012-00
Root Elevators
Fig.12,




Cryer
-53-027-00
Root Elevators
Fig.27,



Cryer 28
-53-028-00
Root Elevators
Fig.28,



Seldin 34
XU-53-034-00
Root Elevators
Fig.34,



Seldin 34S
XU-53-034-09
Root Elevators
Fig.34 S,




Cryer 27S
XU-53-039-00
Root Elevators
Fig. 39,

Root Elevators




Cryer 28S
XU-53-040-00
Root Elevators
Fig. 40,



Flohr
-53-046-00
Root Elevators
Fig. 46,



Barry 50
-53-050-00
Root Elevators
Fig. 50,



Barry 51
-53-051-00
Root Elevators
Fig. 51,



Lindo Levian 1
XU-53-200-01
Root Elevators
Fig. 1,




Lindo Levian 3
XU-53-200-03
Fig. 3,


Root Elevators




Howard-Apexo # 1
XU-53-204-00
Root Elevators
Fig. 1,



XU-53-204-01
Root Elevators
Fig. 1,




Root Elevators
XU-53-205-00
Fig. 1,



Root Elevators
XU-53-205-01
Fig. 1,




Bein 1
-53-220-01
Root Elevator
fig. 2,




Bein 2
-53-220-02
Root Elevator
Fig. 3,


Root Elevators



Bein 3
-53-220-03
Root Elevator
Fig. 3,




Bein 4
-53-220-04
Fig. 4,




Bein 1A
XU-53-221-
01
Fig. 1 A,



Bein 2A
XU-53-221-
02
Fig. 2 A,



Bein 3A
XU-53-221-
03
Fig. 3 A,




Bein 4A
XU-53-221-
04
Fig. 4 A,

Root Elevators




Apical 1
XU-53-222-
01

Fig. 1,




Apical 2
XU-53-222-
02

Fig. 2,




Apical 3
XU-53-222-
03

Fig. 3,




Apical 4
XU-53-222-
04

Fig. 4,



Apical 1A
XU-53-223-
01


Fig. 1 A,



Apical 2A
XU-53-223-
02

Fig. 2 A,

Root Elevators



Apical 3A
XU-53-223-
03

Fig. 3 A,



Apical 4A
XU-53-223-
04

Fig. 4 A,




Apical 301
XU-53-301-00
Root Elevators
Fig. 301,



Apical 302
XU-53-302-00
Root Elevators
Fig. 302,



Apical-Flohr
XU-53-303-00
Root Elevators
Fig. 303,



Apical 304
XU-53-304-00
Root Elevators
Fig. 304,

Root Elevators



Apical 46R
XU-53-311-00
Root Elevators
Fig. 46 R,



Coleman
XU-53-405-01
Root Elevators
Fig. 1,




Coleman 2L
XU-53-405-02
Root Elevators
Fig. 2 L,



Coleman 3R
XU-53-405-03
Root Elevators
Fig. 3 R,




Coupland Bone chisels 1
XU-53-415-01
Root Elevators
Fig. 1,



Coupland Bone Chisels 2
XU-53-415-02
Root Elevators
Fig. 2,

Root Elevators



Coupland Bone chisels 3

XU-53-415-03

Root Elevators
Fig. 3,



Friedman 1

XU-53-505-01

Root Elevators
Fig. 1,



Friedman 2L

XU-53-505-02

Root Elevators
Fig. 2 L,



Friedman 2R

XU-53-505-03


Root Elevators
Fig. 3 R,



Hyllin 1

-53-508-01

Root Elevators
Fig. 1,




Hyllin 2


-53-508-02

Root Elevators
Fig. 2,

Root Elevators




Hylin 3L
-53-508-03
Root Elevators
Fig. 3 L,




Hylin 4R
XU-53-508-04
Root Elevators
Fig. 4 R,



Bein 0
-53-511-00
Root Elevators
Fig. 0,



Bein 1
-53-511-01
Root Elevators
Fig. 1,




Bein 2
-53-511-02
Root Elevators
Fig. 2,




Bein
-53-512-00
Root Elevators
Fig. 0,

Root Elevators




Bein 1
-53-512-01
Root Elevators
Fig. 1,




Bein 2
-53-512-02
Root Elevators
Fig. 2,



Bein 0
-53-513-00
Root Elevators
Fig. 0,



Bein 1
-53-513-01
Root Elevators
Fig. 1,




Bein 2
-53-513-02
Root Elevators
Fig. 2,




Bein 0A
-53-514-00
Root Elevators
Fig. 0 A,

Root Elevators




Bein 1A
-53-514-01
Root Elevators
Fig. 1 A,




Bein 2A
-53-514-02
Root Elevators
Fig. 2 A,



Bein 0A
-53-515-00
Root Elevators
Fig. 0 A,



Bein 1A
-53-515-01
Root Elevators
Fig. 1 A,



Bein 2A
-53-515-02
Root Elevators
Fig. 2 A,



W.F. Barry 320L
XU-53-533-20
Root Elevators
Fig. 320 L,

Root Elevators



W.F. Barry 321R
XU-53-533-21
Root Elevators
Fig. 321 R,



Barry 322
XU-53-533-22
Root Elevators
Fig. 322,



Barry 323
XU-53-533-23
Root Elevators
Fig. 323,



W.F. Barry 370R
XU-53-533-70
Root Elevators
Fig. 370 R,



W.F. Barry 371R
XU-53-533-71
Root Elevators
Fig. 371 R,



Winter 11L
XU-53-551-05
Root Elevators
Fig. 11 L,

Root Elevators



Winter 11R
XU-53-551-08
Root Elevators
Fig. 11 R,



Winter 13L
XU-53-553-05
Root Elevators
Fig. 13 L,



Winter 13R
XU-53-553-08
Root Elevators
Fig. 13 R,



Winter 14L
XU-53-554-05
Root Elevators
Fig. 14 L,



Winter 14R
XU-53-554-08
Root Elevators
Fig. 14 R,



Lecluse 1
XU-53-561-01
Root Elevators
Fig. 1,

Root Elevators



Pott 1
-53-567-01
Root Elevators
Fig. 1,



Pott 2
-53-567-02
Root Elevators
Fig. 2,



Krallenheber 1L
XU-53-571-01
Root Elevators
Fig. 1 L,



Krallenheber 2R
XU-53-571-02
Root Elevators
Fig. 2 R,



Warwick-James 1
XU-53-576-01
Root Elevators
Fig. 1,



Warwick-James 2L
XU-53-576-02
Root Elevators
Fig. 2 L,

Root Elevators




Warwick-James 3R
XU-53-576-03
Root Elevators
Fig. 3 R,



Chompret 1
XU-53-577-01
Root Elevators
Fig. 1,



Chompret 2
XU-53-577-02
Root Elevators
Fig. 2,



Chompret 3
XU-53-577-03
Root Elevators
Fig. 3,




Heidebrink 1
XU-53-582-01
Root Elevators
Fig. 1,



Heidebrink 2
XU-53-582-02
Root Elevators
Fig. 2,


Root Elevators




Heidebrink 3
XU-53-582-03
Root Elevators
Fig. 3,



West 1
-53-584-01
Root Elevators
Fig. 1,



West 2
-53-584-02
Root Elevators
Fig. 2,



West 3
-53-584-03
Root Elevators
Fig. 3,

Bone Rongeurs



Bone Rongeur

XU-54-101

Bone Holding Forceps
185 mm, 7 1/2",



Bone Rongeur

XU-54-102

Bone Holding Forceps
185 mm, 7 1/2",



Bone Rongeur

XU-54-103


Bone Holding Forceps
200 mm, 8",



Bone Rongeur

XU-54-104

Bone Holding Forceps
215 mm, 8 1/2",



Bone Rongeur

XU-54-105

Bone Rongeurs
145 mm, 5 3/4",



Friedman-Micro 6"

XU-54-106

Bone Rongeurs
150 mm, 6",

Bone Rongeurs



Friedman-Micro 5"

XU-54-107

Bone Rongeurs
145 mm, 5\|\|, \|\|



Bone Rongeur

XU-54-108

Bone Rongeurs
140 mm, 5?", \|\|



LUER 6"

-54-109

Bone Rongeurs
150 mm, 6,



LUER 6"

-54-110

Bone Rongeurs
155 mm, 6,



2A-MEAD, SQUARE TIPS 5"

XU-54-111

Bone Rongeur
145 mm, 5,



BLUMENTHAL 6"

XU-54-112

Bone Rongeur
155 mm, 6,

Bone Rongeurs



BLUMENTHAL 45° 6"

XU-54-113

Bone Rongeur
150 mm, 6"



Bone Rongeur 7"

XU-54-114

Bone Rongeur
175 mm, 7"



Bone Rongeur

XU-54-115

Bone Rongeur
170 mm, 6 7/8"



BEYER 7"

XU-54-116

Bone Rongeur
180 mm, 7"



MEAD 6 1/2"

XU-54-117

Bone Rongeur
165 mm, 6 1/2"



BANE-HEARTMANN 7"

XU-54-118

Bone Rongeur
175 mm, 7"

Bone Rongeurs



BEYER DOUBLE-ACTION 7"

XU-54-119

Bone Rongeur
180 mm, 7",



ZUFAL-JANSEN 7"

XU-54-120

Bone Rongeur
180 mm, 7",



Bone Rongeur

XU-54-121

Bone Rongeur
130 mm, 5½",



Bone Rongeur

XU-54-122

Bone Rongeur
130 mm, 5½",



Cleveland 5"

XU-54-123

Bone Cutting Forceps
130 mm, 5",



Bone Rongeur

XU-54-124

Bone Cutting Forceps
150 mm, 6",

Bone Files



Fig. 1
Polokoff

Polokoff
XU-55-101
Bone Files
Fig- 1,




Fig. 33

XU-55-102
Bone Files
Fig- 33,




Fig. 64
Miler

Miler
XU-55-103
Bone Files
Fig- 64,



Fig. 45
Miler

Miler
XU-55-104
Bone Files
Fig- 45,



Fig. 10

XU-55-105
Bone Files
Fig- 10,



Fig. 1
Wahl

Wahl
XU-55-106
Bone Files
Fig- 1,

Bone Files




Wahl
XU-55-107
Bone Files
Fig- 2,



XU-55-108
Bone Files
Fig- 12CA,



XU-55-109
Bone Files
Fig- 12A,



Howard
XU-55-110
Bone Files
Fig- 12,

Bone Currettes



Fig. 1
Hemingway

Hemingway
XU-56-101
Bone Currettes
Fig- 1,



Fig. 2
Hemingway

Hemingway
XU-56-102
Bone Currettes
Fig- 2,



Fig. 3
Hemingway

Hemingway
XU-56-103
Bone Currettes
Fig- 3,



Fig. 4
Hemingway

Hemingway
XU-56-104
Bone Currettes
Fig- 4,



Fig. 85
Lucas

Lucas
-56-105
Bone Currettes
Fig- 85,



Fig. 86
Lucas

Lucas
-56-106
Bone Currettes
Fig- 86,

Bone Curettes



Fig. 87
Lucas

Lucas
-56-107
Bone Curettes
Fig- 87,



Fig. 88
Lucas

Lucas
-56-108
Bone Curettes
Fig- 88,



Fig. 1
Hemingway

Hemingway
XU-56-109
Bone Curettes
Fig- 1,



Fig. 2
Hemingway

Hemingway
XU-56-110
Bone Curettes
Fig- 2,



Fig. 3
Hemingway

Hemingway
XU-56-111
Bone Curettes
Fig- 3,



Fig. 4
Hemingway

Hemingway
XU-56-112
Bone Curettes
Fig- 4,

Bone Currettes



Fig. 85
Lucas

Lucas

-56-113

Bone Currettes
Fig- 85,



Fig. 86
Lucas

Lucas

-56-114

Bone Currettes
Fig- 86,



Fig. 87
Lucas

Lucas

-56-115

Bone Currettes
Fig- 87,



Fig. 88
Lucas

Lucas

-56-116

Bone Currettes
Fig- 88,

Cotton & Dressing Tweezers



Cotton & Dressing Tweezers

XU-57-101

Tooth Tweezers
130 mm, 5 $\frac{1}{8}$ ",



College Tweezer

XU-57-102

Cotton & Dressing Tweezers
130 mm, 5 $\frac{1}{8}$ ",



College Tweezer Smooth Plan Handle

XU-57-103

Cotton & Dressing Tweezers
150 mm, 6",



Meriam Tweezers

XU-57-104

Cotton & Dressing Tweezers
150 mm, 6",



College Tweezer

XU-57-105

Cotton & Dressing Tweezers
160 mm, 6",



Meriam Tweezers

XU-57-106

Cotton & Dressing Tweezers
160 mm, 6",

Cotton & Dressing Tweezers



Meriam Tweezers

XU-57-107

Cotton & Dressing Tweezers
160 mm, 6,



College Tweezer

XU-57-108

Cotton & Dressing Tweezers
150 mm, 6,



College Tweezer

XU-57-109

Cotton & Dressing Tweezers
150 mm, 6,



Self Locking Tweezer Serr.

XU-57-110

Cotton & Dressing Tweezers
150 mm, 6,



College Tweezer

XU-57-111

Cotton & Dressing Tweezers
150 mm, 6,



College Tweezer

XU-57-112

Cotton & Dressing Tweezers
150 mm, 6,

Cotton & Dressing Tweezers



Perry Tweezer

XU-57-113

Cotton & Dressing Tweezers
130 mm,



XU-57-114

Cotton & Dressing Tweezers
150 mm,6,



Flagg Tweezer Smooth

XU-57-115

Cotton & Dressing Tweezers
150 mm,6,



Flagg Tweezer Serrated

XU-57-116

Cotton & Dressing Tweezers
150 mm,6,



XU-57-117

Cotton & Dressing Tweezers
150 mm,6,



Cotton & Dressing Tweezers

XU-57-118

Tooth Tweezers
150 mm,6",



Cotton & Dressing Tweezers

XU-57-119

Tooth Tweezers
150 mm, 6",



Self Locking Tweezer Smooth

XU-57-120

Cotton & Dressing Tweezers
150 mm, 6",



Suture Forcep

XU-57-121

Cotton & Dressing Tweezers
150 mm, 6",

Rubber Dam Punch & Clamps



Rubber Dam Punch Ivory
XU-58-101
Rubber-Dam Instruments
160 mm, 6,



Anisworth
XU-58-102
Rubber-Dam Instruments
160 mm, 6?,



Brewer
-58-103
Rubber-Dam Instruments
170 mm, 6?,



Palmer
-58-104
Rubber-Dam Instruments
175 mm, 7,

Orthodontic Pliers



BIRD BEAK PLIER
XU-59-101
ORTHODONTIC PLIERS, ANGLE



CUTTING PLIER FOR ORTHODONTICS
XU-59-102
ORTHODONTIC PLIERS NO.1



CUTTING PLIERS FOR ORTHODONTIC
XU-59-103
ORTHODONTIC PLIERS NO.2



DE LA ROSA ORTHO PLIER
XU-59-104
ORTHODONTIC GROOVED PLIERS, DE LA
ROSA



DE LA ROSA ORTHO PLIERS
XU-59-105
ORTHODONTIC NON GROOVED PLIERS,
DE LA ROSA



DISTAL END CUTTER
XU-59-106
DISTAL END CUTTER

Orthodontic Pliers



DISTAL END CUTTER
XU-59-107 ITC



HARD WIRE CUTTER TC
XU-59-108 3TC



LIGHT WIRE CUTTER
XU-59-109
MICRO-MINIATURE WIRE CUTTER



LIGHT WIRE CUTTER
XU-59-110 ITC



MC KELLOPS PLIERS FOR ORTHODON
XU-59-111
ORTHODONTIC PLIERS, MC KELLOPS



MICRO DISTAL END CUTTER
XU-59-112
MICRO DISTAL WIRE CUTTER

Orthodontic Pliers



NANCE PLIERS FOR ORTHODONTICS

XU-59-113

ORTHODONTIC PLIERS NANCE



ORTHODONTIC PLIERS

XU-59-114

ORTHODONTIC PLIERS TO FORM
LINGUAL ARCHS



ORTHODONTIC PLIER

XU-59-115

ORTHODONTIC PLIERS TO REMOVE
BRACKETS



ORTHODONTIC PLIER

XU-59-116

ORTHODONTIC ATTACHMENT REMOVER
PLIERS



ORTHODONTIC PLIERS

XU-59-117

ORTHODONTIC ATTACHMENTS
REMOVER PLIERS



ORTHODONTIC PLIERS

XU-59-118

ORTHODONTIC DIRECT BONDING PLIERS

Orthodontic Pliers



PLIERS FOR ORTHODONTIC CURVED
XU-59-119
ORTHODONTIC PLIERS HOW NO.2
CURVED



PLIERS FOR ORTHODONTICS
XU-59-120
ORTHODONTIC PLIERS



PLIERS FOR ORTHODONTICS
XU-59-121
ORTHODONTIC PLIERS NO.2



PLIERS FOR ORTHODONTICS
XU-59-122
ORTHODONTIC PLIERS WITH SMOOTH
BEAKS



PLIERS FOR ORTHODONTICS
XU-59-123
ORTHODONTIC PLIERS NO.2



PLIERS FOR ORTHODONTICS
XU-59-124
ORTHODONTIC CONTOURING PLIERS,
REYNOLDS

Orthodontic Pliers



PLIERS FOR ORTHODONTICS
XU-59-125
ORTHODONTIC PLIERS, WALDSACHS



PLIERS FOR ORTHODONTICS
XU-59-126
ORTHODONTIC PLIERS JOHNSON



PLIERS FOR ORTHODONTICS
XU-59-127
ORTHODONTIC PLIERS, NANCE



PLIERS FOR ORTHODONTICS
XU-59-128
ORTHODONTIC PLIERS, YOUNG



PLIERS FOR ORTHODONTICS
XU-59-129
ORTHODONTIC PLIERS



PLIERS FOR ORTHODONTICS
XU-59-130
ORTHODONTIC PLIERS WITH ROUNDED
JAW

Orthodontic Pliers



PLIERS FOR ORTHODONTICS
XU-59-131
ORTHODONTIC PLIERS, ANGLE



PLIERS FOR ORTHODONTICS
XU-59-132
ORTHODONTIC PLIERS, ADAMS



PLIERS FOR ORTHODONTICS
XU-59-133
ORTHODONTIC PLIERS WITH ROUND
NOSE



PLIERS FOR ORTHODONTICS
XU-59-134
ORTHODONTIC PLIERS, ADERER



PLIERS FOR ORTHODONTICS
XU-59-135
ORTHODONTIC PLIERS ADERER



PLIERS FOR ORTHODONTICS
XU-59-136
ORTHODONTIC PLIERS, TWEED

Orthodontic Pliers



PLIERS FOR ORTHODONTICS

XU-59-137

ORTHODONTIC PLIERS



PLIERS FOR ORTHODONTICS

XU-59-138

ORTHODONTIC PLIERS, ANGLE NO.1



PLIERS FOR ORTHODONTICS STR.

XU-59-139

ORTHODONTIC PLIERS HOW NO.1
STRAIGHT



POSTERIOR BAND REMOVER

XU-59-140

ORTHODONTIC BAND REMOVER PLIERS



POSTERIOR BAND REMOVER

XU-59-141

ORTHODONTIC BAND REMOVER PLIERS



UNIVERSAL PLIERS FOR ORTHODON.

XU-59-142

UNIVERSAL ORTHODONTIC PLIERS

Orthodontic Pliers



UNIVERSAL ORTHO PLIERS

XU-59-143

UNIVERSAL BENDING ORTHODONTIC
PLIERS



ROUND AND CONCAVE PLIERS

XU-59-144

ORTHODONTIC PLIERS, ROUND AND
CONCAVE



TWEED ORTHO PLIER

XU-59-145

ORTHODONTIC PLIERS, TWEED



TWEED ORTHO PLIERS

XU-59-146

ORTHODONTIC PLIERS, TWEED



TWEED ORTHO PLIERS

XU-59-147

ORTHODONTIC LOOP FORMING
PLIERS, TWEED



WEINGART ORTHO PLIERS

XU-59-148

ORTHODONTIC UNIVERSAL PLIERS,
WEINGART

Orthodontic Pliers



WIRE CUTTER

XU-59-151

WIRE CUTTER FOR ANY SIZE LIGATURE
WIRES



WIRE CUTTER

XU-59-152 ITC



HARD WIRE CUTTER TC

XU-59-153-3TC

Needle Holders



Crile-Murray
XU-60-101
Needle Holders
130 mm, 5,



Mayo-Hegar
XU-60-102
Needle Holders
145 mm, 5,



Crile-Murray
XU-60-103
Needle Holders
150 mm, 6,



Mayo-Hegar
XU-60-104
Needle Holders
150 mm, 6,




Mayo-Hegar
XU-60-105
Needle Holders
180 mm, 7,



Mayo-Hegar
XU-60-106
Needle Holders
200 mm, 8,

Needle Holders



Olsen-Hegar
XU-60-107
Needle Holders
140 mm, 5,



Gillies
-60-108
Needle Holders
150 mm, 6,



Boynton
-60-109
Needle Holders
125 mm, 5,



Mathieu
-60-110
Needle Holders
170 mm, 6,



Mathieu
-60-111
Needle Holders
200 mm, 8,



Crile
-60-112
Needle Holders
150 mm, 6,

Needle Holders



Toennis
-60-113
Needle Holders
180 mm, 7,



Mayo-Hegar
XU-60-114
Needle Holders
165 mm, 6,



Mayo-Hegar
XU-60-115
Needle Holders
210 mm, 8,



Mayo-Hegar
XU-60-116
Needle Holders
165 mm, 6,



Mayo-Hegar
XU-60-117
Needle Holders
210 mm, 8,



Debaquey
-60-118
Needle Holders
150 mm, 6,

Needle Holders



Debakey
-60-119
Needle Holders
180 mm, 7,



Ryder
-60-120
Needle Holders
135 mm, 5,



Ryder
-75-121
Needle Holders
155 mm, 6,



Ryder 7 CARB-N-SERT Jaws
XU-60-122
Needle Holders
175 mm, 7,



Bozemann 5
XU-60-123
Needle Holders
135 mm, 5,



Olsen-Hegar CARB-N-SERT Jaws
XU-60-124
Needle Holders
125 mm, 5,

Needle Holders



Olsen-Hegar CARB-N-SERT Jaws

XU-60-125

Needle Holders
165 mm, 6,



Olsen-Hegar CARB-N-SERT Jaws

XU-60-126

Needle Holders
160 mm, 6,



Olsen-Hegar CARB-N-SERT Jaws

XU-60-127

Needle Holders
145 mm, 5,



Needle Holders

XU-60-128

Needle Holders
185 mm, 7 1/2",



Needle Holders

XU-60-129

Needle Holders
180 mm, 7",



Ryder 8 CARB-N-SERT Jaws

XU-60-130

Needle Holders
200 mm, 8,

Needle Holders



Needle Holders
XU-60-131
Needle Holders
200 mm, 8",



Mathieu 5
XU-60-132
Needle Holders
140 mm, 5,



Mathieu 6
XU-60-133
Needle Holders
170 mm, 6,



Mathieu 8
XU-60-134
Needle Holders
200 mm, 8,



Mathieu 5
XU-60-135
Needle Holders
140 mm, 5,



Mathieu 6
XU-60-136
Needle Holders
170 mm, 6,

Needle Holders



Mathieu 8
XU-60-137
Needle Holders
200 mm, 8,



Mathieu 7
XU-60-138
Needle Holders
195 mm, 7,



Mathieu CARB-N-SERT Jaws
XU-60-139
Needle Holders
205 mm, 8,

Gum Scissors



LISTER Bandage Scissors

XU-61-101

GUM Scissors
Extra fine
4-1/2" (11.4 cm),



LISTER Bandage Scissors

XU-61-102

Gum Scissors
Extra fine
5-1/2" (14 cm),



LISTER Bandage Scissors

XU-61-103

Gum Scissors
Extra fine
7-1/4" (18.4 cm),



KNOWLES Bandage Scissors

XU-61-104

Gum Scissors
Straight
5-1/2" (14 cm),



HEATH Suture Scissors

XU-61-105

Gum Scissors
Serrated blade
6-1/8" (15.6 cm),



SPENCER Stitch Scissors

XU-61-106

Gum Scissors
Delicate
3-1/2" (8.9 cm),

Gum Scissors



SHORTBENT Stitch Scissors

XU-61-107

Gum Scissors
Curved delicate
3-1/2" (8.9 cm),



LITTAUER JR. Stitch Scissors

XU-61-108

Gum Scissors
4-1/2" (11.4 cm),



NORTHBENT Stitch Scissors

XU-61-109

Gum Scissors
Curved light pattern
4-3/4" (12.1 cm),



LITTAUER Stitch Scissors

XU-61-110

Gum Scissors
Standard pattern
5-1/2" (14 cm),



LONG Oral Surgery Stitch Scissors

XU-61-111

Gum Scissors
6" (15.2 cm),



MILTEX Stitch Scissors

XU-61-112

Gum Scissors
Angled delicate
4-1/2" (11.4 cm),

Gum Scissors



REEH Stitch Scissors

XU-61-113

Gum Scissors

Sharp tips small hook on 1 blade ribbon type
3-3/4" (9.5 cm),



MILTEX Stitch Scissors

XU-61-114

Gum Scissors

Curved sharp pointed tips ribbon type
3-7/8" (9.8 cm),



MILTEX Stitch Scissors

XU-61-115

Gum Scissors

Curved sharp pointed tips ribbon type
4-1/2" (11.4 cm),



VANTAGE Crown Scissors

XU-61-116

Gum Scissors

Collar curved
4" (10.2 cm),



VANTAGE Collar

XU-61-117

Gum Scissors

Crown Scissors straight
4" (10.2 cm),



MILTEX Wire Cutting Scissors

XU-61-118

Gum Scissors

Angled to side one serrated blade
4-3/4" (12.1 cm),

Gum Scissors



MILTEX Wire Cutting Scissors

XU-61-119

Gum Scissors

Angled to side, one serrated blade, Carb-N-Sert inserts
4-3/4" (12.1 cm),



Standard Pattern Operating Scissors

XU-61-120

Gum Scissors

Straight sharp-sharp points
4-1/2" (11.4 cm),



Standard Pattern Operating Scissors

XU-61-121

Gum Scissors

Straight sharp-sharp points
5-1/2" (14 cm),



Standard Pattern Operating Scissors

XU-61-122

Gum Scissors

Straight sharp-blunt points
4-1/2" (11.4 cm),



Standard Pattern Operating Scissors

XU-61-123

Gum Scissors

Straight sharp-blunt points
5-1/2" (14 cm),



Standard Pattern Operating Scissors

XU-61-124

Gum Scissors

Straight sharp-blunt points
5-1/2" (14 cm),

Gum Scissors



Standard Pattern Operating Scissors

XU-61-125

Gum Scissors
Curved sharp/sharp points
4-1/2" (11.4 cm),



Standard Pattern Operating Scissors

XU-61-126

Gum Scissors
Curved sharp/sharp points
5-1/2" (14 cm),



Standard Pattern Operating Scissors

XU-61-127

Gum Scissors
Curved
4-1/2,



Standard Pattern Operating Scissors

XU-61-128

Gum Scissors
Curved, sharp-blunt points
5-1/2" (14 cm),



MAYO Dissecting Scissors

XU-61-129

Gum Scissors
Straight, standard beveled blades
5-1/2" (14 cm),



MAYO Dissecting Scissors

XU-61-130

Gum Scissors
Straight, standard beveled blades
5-1/2" (14 cm),

Gum Scissors



MAYO Dissecting Scissors

XU-61-131

Gum Scissors
Curved, standard beveled blades
5-1/2" (14 cm),



MAYO Dissecting Scissors

XU-61-132

Gum Scissors
Curved, standard beveled blades
5-1/2" (14 cm),



Delicate Pattern METZENBAUM Scissors

XU-61-133

Gum Scissors
Straight, blunt points
5-1/2" (14 cm),



Delicate Pattern METZENBAUM Scissors

XU-61-134

Gum Scissors
Straight, sharp points
5-1/2" (14 cm),



Standard Pattern METZENBAUM Scissors

XU-61-135

Gum Scissors
Curved, blunt points
5-1/2" (14 cm),



Standard Pattern METZENBAUM Scissors

XU-61-136

Gum Scissors
Curved, blunt points
5-1/2" (14 cm),

Gum Scissors



Delicate Pattern METZENBAUM Scissors

XU-61-137

Gum Scissors
Curved, blunt points
5-1/2" (14 cm),



Delicate Pattern METZENBAUM Scissors

XU-61-138

Gum Scissors
Curved, sharp points
5-1/2" (14 cm),



KELLY Scissors

XU-61-139

Gum Scissors
Straight, sharp points
6-1/4" (15.9 cm),



KELLY Scissors

XU-61-140

Gum Scissors
Curved, sharp points
6-1/4" (15.9 cm),



KELLY Scissors

XU-61-141

Gum Scissors
Sharp points, angled on side
6-1/4" (15.9 cm),



KAYE Blepharoplasty and Dissection Scissors

XU-61-142

Gum Scissors
One serrated blade, ribbon type
4-1/4" (10.8 cm),

Gum Scissors



DEAN Dissecting Scissors

XU-61-143

Gum Scissors

Blades angled on flat, one serrated blade
6-3/4" (17.1 cm),



KAYE Dissecting Scissors

XU-61-144

Gum Scissors

Curved, one serrated blade
4-17/32" (11.5 cm),



WAGNER Plastic Surgery Scissors

XU-61-145

Gum Scissors

Straight, sharp-sharp points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-146

Gum Scissors

Straight, sharp-sharp points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-147

Gum Scissors

Straight, sharp-blunt points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-148

Gum Scissors

Straight, sharp-blunt points, serrated blade
4-3/4" (12.1 cm),

Gum Scissors



WAGNER Plastic Surgery Scissors

XU-61-149

Gum Scissors
Straight, blunt-blunt points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-150

Gum Scissors
Straight, blunt-blunt points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-151

Gum Scissors
Curved, sharp-sharp points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-152

Gum Scissors
Curved, sharp-sharp points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-153

Gum Scissors
Curved, sharp-blunt points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-154

Gum Scissors
Curved, sharp-blunt points, serrated blade
4-3/4" (12.1 cm),

Gum Scissors



WAGNER Plastic Surgery Scissors

XU-61-155

Gum Scissors
Curved, blunt-blunt points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-156

Gum Scissors
Curved, blunt-blunt points, serrated blade
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-157

Gum Scissors
Straight, sharp-sharp points
4-3/4" (12.1 cm),



WAGNER Plastic Surgery Scissors

XU-61-158

Gum Scissors
Straight, sharp-sharp point, serrated blade
4-3/4" (12.1 cm),



Delicate Pattern METZENBAUM Scissors

XU-61-159

Gum Scissors
Curved
5" (12.7 cm),



Iris Scissors

XU-61-160

Gum Scissors
Straight
4" (10.2 cm),

Gum Scissors



Iris Scissors
XU-61-161
Gum Scissors
Curved
4" (10.2 cm),



Iris Scissors
XU-61-162
Gum Scissors
Straight
4-1/2" (11.4 cm),



Iris Scissors
XU-61-163
Gum Scissors
Straight
4-1/2" (11.4 cm),



Iris Scissors
XU-61-164
Gum Scissors
Curved
4-1/2" (11.4 cm),



Iris Scissors
XU-61-165
Gum Scissors
Curved
4-1/2" (11.4 cm),



Iris Scissors
XU-61-166
Gum Scissors
Angled on side
4-1/2" (11.4 cm),

Gum Scissors



BROPHY Scissors

XU-61-167

Gum Scissors
Straight, sharp points
5-1/2" (14 cm),



BROPHY Scissors

XU-61-168

Gum Scissors
Curved, sharp points
5-1/2" (14 cm),



GOLDMAN-FOX Scissors

XU-61-169

Gum Scissors
Curved fine tips, one serrated blade
4" (10.2 cm),



GOLDMAN-FOX Scissors

XU-61-170

Gum Scissors
Curved fine tips, one serrated blade
5-1/2" (14 cm),



KELLY Scissors

XU-61-171

Gum Scissors
Sharp points, straight, one serrated blade
6-1/4" (15.9 cm),



KELLY Scissors

XU-61-172

Gum Scissors
Sharp points, curved, one serrated blade
6-1/4" (15.9 cm),

Gum Scissors



KELLY Scissors

XU-61-173

Gum Scissors

Sharp points, angular, one serrated blade
6-1/4" (15.9 cm),



Iris Operating Scissors

XU-61-174

Gum Scissors

Curved sideways
4-1/2" (11.4 cm),



GOLDMAN-FOX Scissors

XU-61-175

Gum Scissors

Curved, one serrated blade, angled shank
4-1/2" (11.4 cm),



LA GRANGE Scissors

XU-61-176

Gum Scissors

Curved, one serrated blade, angled shank
4-1/4" (10.8 cm),



BROPHY Operating Scissors

XU-61-177

Gum Scissors

Curved, sharp/sharp, with one serrated blade
5-1/2" (14 cm),



LOCKLIN Operating Scissors

XU-61-178

Gum Scissors

Angled on side 25 degrees, one serrate
blade
6-1/4" (15.9 cm),

Gum Scissors



LOCKLIN Operating Scissors

XU-61-179

Gum Scissors

Curved shank, angled on side 25 degrees,
one serrated blade
6-1/4" (15.9 cm),



MILTEX SuperCut GOLDMAN-FOX Gum Scissors

XU-61-180

Gum Scissors

Curved
4-1/2" (11.4 cm),



HEATH Suture Scissors

XU-61-181

Gum Scissors

Serrated blade
6-1/8" (15.6 cm),



Iris Scissors

XU-61-182

Gum Scissors

Straight, with 20 mm blades, delicate
3-1/2" (8.9 cm),



Iris Scissors

XU-61-183

Gum Scissors

Curved, with 20 mm blades, delicate
3-1/2" (8.9 cm),



NOYES (ROSS) Iris Scissors

XU-61-184

Gum Scissors

Curved, sharp/sharp points
4-1/2" (11.4 cm),

Gum Scissors



LOCKLIN Operating Scissors

XU-61-179

Gum Scissors

Curved shank, angled on side 25 degrees,
one serrated blade
6-1/4" (15.9 cm),



MILTEX SuperCut GOLDMAN-FOX Gum Scissors

XU-61-180

Gum Scissors

Curved
4-1/2" (11.4 cm),



HEATH Suture Scissors

XU-61-181

Gum Scissors

Serrated blade
6-1/8" (15.6 cm),



Iris Scissors

XU-61-182

Gum Scissors

Straight, with 20 mm blades, delicate
3-1/2" (8.9 cm),



Iris Scissors

XU-61-183

Gum Scissors

Curved, with 20 mm blades, delicate
3-1/2" (8.9 cm),



NOYES (ROSS) Iris Scissors

XU-61-184

Gum Scissors

Curved, sharp/sharp points
4-1/2" (11.4 cm),

Gum Scissors



CASTROVIEJO Corneal Scissors

XU-61-185

Gum Scissors
Curved, sharp points
3-3/4" (9.5 cm),



QUIMBY Gum Scissors

XU-61-186

Gum Scissors
Strongly curved and pointed, smooth blade
15 mm
4" (10.2 cm),



JOSEPH Scissors

XU-61-187

Gum Scissors
Straight
5-3/4" (14.6 cm),



JOSEPH Scissors

XU-61-188

Gum Scissors
Curved
5-3/4" (14.6 cm),

Amalgam Instruments



Fig. SHO-A

XU-62-117

Amalgam Instruments
Fig- SHO-A,



Fig. 3-SHO

XU-62-118

Amalgam Instruments
Fig- 3-SHO,



Fig. 104

XU-62-119

Amalgam Instruments
Fig- 104,



Fig. 8

XU-62-120

Amalgam Instruments
Fig- 8,



Fig. 1

XU-62-121

Amalgam Instruments
Fig- 1,



Fig. 2

XU-62-122

Amalgam Instruments
Fig- 2,

Amalgam Instruments



Fig 3

XU-62-123

Amalgam Instruments
Fig- 3,



Fig 4

XU-62-124

Amalgam Instruments
Fig- 4,



Fig. SHO-A

XU-62-125

Amalgam Instruments
Fig- SHO-A,



Fig. 3-SHO

XU-62-126

Amalgam Instruments
Fig- 3-SHO,



Fig. 104

XU-62-127

Amalgam Instruments
Fig- 104,



Fig. 8

XU-62-128

Amalgam Instruments
Fig- 8,

Amalgam Instruments



XU-62-129

Amalgam Instruments
Fig- 1,



XU-62-130

Amalgam Instruments
Fig- 2,



XU-62-131

Amalgam Instruments
Fig- 3,



XU-62-132

Amalgam Instruments
Fig- 4,



XU-62-133

Amalgam Instruments
Fig- 1,



XU-62-134

Amalgam Instruments
Fig- 2,

Amalgam Instruments



Fig. 3

XU-62-135

Amalgam Instruments
Fig- 3,



Fig. 4

XU-62-136


Amalgam Instruments
Fig- 4,



Tofflemire-Universal

XU-62-137

Amalgam Instruments
Ivory Fig- 8,5 mm,6 mm,7 mm,



Nyström I + II

Ivory Fig. 1

Ivory Fig. 12

XU-62-142

Amalgam Instruments
Ivory Fig- 1,Ivory Fig- 2,



Amalgam Instruments

XU-62-143

Amalgam Pistolen
210 mm,



Amalgam Instruments

XU-62-144

Amalgam Pistolen
185 mm,

Amalgam Instruments



Amalgam Instruments

XU-62-145

Amalgam Pistolen
135 mm,



Amalgam Instruments

XU-62-146

Amalgam Pistolen
145 mm,



Amalgam Carrier

XU-62-147

?
Fig. 1, Fig. 2, Fig. 3,



Amalgam Carrier

XU-62-148

?
Fig. 1, Fig. 2, Fig. 3,



Amalgam Carrier

XU-62-149

?
Fig. 1, Fig. 2, Fig. 3,

Mouth Gags & Retractor



Heister
XU-63-141
Mouth Gags,
130 mm, 5 1/8,



MOLT PEDIATRIC
XU-63-142
Mouth Gags,
100 mm, 4,



MOLT CHILD
XU-63-143
Mouth Gags,
120 mm, 4,



MOLT ADULT
XU-63-144
Mouth Gags,
140 mm, 5,



FERGUSON-ACKLAND
XU-63-145
Mouth Gags,
145 mm, 5,



DENHART
XU-63-146
Mouth Gags,
125 mm, 5,

Mouth Gags & Retractor



ROSER-KONIG

XU-63-147

Mouth Gags,
160 mm, 6 1/4,



ROSER-KONIG

XU-63-148

Mouth Gags,
185 mm, 7,



LANGENBECK

XU-63-149

Finger Protectors
185 mm, 7,



XU-63-150

Tongue Depressor and Cheek Retractors



Retractor

XU-63-151

Mouth Gags,



HARTMANN

XU-63-152

Retractors
145 mm, 5,

Mouth Gags & Retractor



Retractor
XU-63-153
Mouth Gags,
140 mm,5,



SCHMID
XU-63-154
Mouth Gags,
140 mm,5,



Retractor
XU-63-155
Mouth Props
140 mm,5?,



Retractor
XU-63-156
Retractors
160 mm,6?,



VESTIBULUM RETRACTOR
XU-63-157
Retractors
140 mm,5,



STERNBERG
XU-63-158
Retractors
120 mm,4,

Mouth Gags & Retractor



ROUX

-63-159

Retractors

19x22 mm, 23x30 mm,



ROUX

-63-160

Retractors

25x26 mm, 28x38 mm,



ROUX

-63-161

Retractors

28x29 mm, 33x45 mm,



MIDDELDORPF

XU-63-162

Retractors

15x15, 200 mm, 8,



MIDDELDORPF

XU-63-163

Retractors

22x22 mm, 210 mm, 8,



MIDDELDORPF

XU-63-164

Cheek Retractors

26x30 mm, 230 mm, 9,

Mouth Gags & Retractor



Gauze Case
XU-63-165
Gauze Case
26x30 mm, 230 mm, 9,



Retractor
XU-63-166
Retractors
165 mm, 6",



Retractor
XU-63-167
Retractors
165 mm, 6",



Retractor
XU-63-168
Retractors
165 mm, 6",



Retractor
XU-63-169
Retractors
165 mm, 6",



Retractor
XU-63-170
Retractors
165 mm, 6",

Mouth Gags & Retractor



Retractor
XU-63-171
Retractors
165 mm, 6 7/8",



Retractor
XU-63-172
Retractors
165 mm, 6 7/8",



Retractor
XU-63-173
Retractors
165 mm, 6 7/8",



Retractor
XU-63-174
Retractors
160 mm, 6 3/8",



Retractor
XU-63-175
Retractors
160 mm, 6 3/8",



LANGENBECK-GREEN
XU-63-176
Retractors
160 mm, 16x6 mm,

Periodontia Instruments



Fig. 1/2
105/100
Kirkland

Kirkland

-64-101

Periodontia Instruments

Fig- 1/2,



Fig. 1/2
Orban

Orban

-64-102

Periodontia Instruments

Fig- 1/2,



Fig. 1/2
105/100
Kirkland

Kirkland

-64-103

Periodontia Instruments

Fig- 1/2,



Fig. 1/2
Orban

Orban

-64-104

Periodontia Instruments

Fig- 1/2,



Fig. 7
Goldman-Fox

Goldman-Fox

XU-64-105

Periodontia Instruments

Fig- 7,



Fig. 8
Goldman-Fox

Goldman-Fox

XU-64-106

Periodontia Instruments

Fig- 8,

Periodontia Instruments



Fig. 9
Goldman-Fox

Goldman-Fox
XU-64-107
Periodontia Instruments
Fig- 9,



Fig. 11
Goldman-Fox

Goldman-Fox
XU-64-108
Periodontia Instruments
Fig- 11,



Fig. 3/4
Buck

Buck
-64-109
Periodontia Instruments
Fig- 3/4,



Fig. 5/6
Buck

Buck
-64-110
Periodontia Instruments
Fig- 5/6,



Fig. 7
Goldman-Fox

Goldman-Fox
XU-64-111
Periodontia Instruments
Fig- 7,



Fig. 8
Goldman-Fox

Goldman-Fox
XU-64-112
Periodontia Instruments
Fig- 8,

Periodontia Instruments



Fig. 9
Goldman-Fox

Goldman-Fox
XU-64-113
Periodontia Instruments
Fig- 9,



Fig. 11
Goldman-Fox

Goldman-Fox
XU-64-114
Periodontia Instruments
Fig- 11,



Fig. 3/4
Buck

Buck
-64-115
Periodontia Instruments
Fig- 3/4,



Fig. 5/6
Buck

Buck
-64-116
Periodontia Instruments
Fig- 5/6,



Fig. 1/2
Sanders

Sanders
-64-117
Periodontia Instruments
Fig- 1/2,



Fig. 4/5
Sanders

Sanders
-64-118
Periodontia Instruments
Fig- 4/5,

Periodontia Instruments



Fig. 19/20

OCHSENBEIN CHISEL NO 4

XU-64-119

Periodontia Instruments
Fig- 19/20,



Fig. 5/6
Back

Crane-Kaplan

XU-64-120

Periodontia Instruments
Fig- CK3,



Fig. 1/2
Back view

Sanders

-64-121

Periodontia Instruments
Fig- 1/2,



Fig. 4/5
Sanders

Sanders

-64-122

Periodontia Instruments
Fig- 4/5,



Fig. 19/20

OCHSENBEIN CHISEL NO 4

XU-64-123

Periodontal Instruments
Fig- 19/20,



Fig. CK3
Graber's radius

Crane-Kaplan

XU-64-124

Periodontia Instruments
Fig- CK3,

Periodontia Instruments



Fig. HU 1/2
Curette

Curette

-64-125

Periodontia Instruments
Fig- HU 1/2,



Fig. 2R/2L
Columbia

Columbia

-64-126

Periodontia Instruments
Fig- 2R/2L,



Fig. 4R/4L
Columbia

Columbia

-64-127

Periodontia Instruments
Fig- 4R/4L,



Fig. 13/14
Columbia

Columbia

-64-128

Periodontia Instruments
Fig- 13/14,



Fig. 411/412

XU-64-129

Periodontia Instruments
Fig- 411/412,



Fig. 413/414

XU-64-130

Periodontia Instruments
Fig- 413/414,

Periodontia Instruments



Zerling-Chiles

Zerling-Chiles

XU-64-131

Periodontia Instruments

Fig- 413/414,



Fig. HU 1/2
Curette

Curette

-64-132

Periodontia Instruments

Fig- HU 1/2,



Fig. 2R/2L
Columbia

Columbia

-64-133

Periodontia Instruments

Fig- 2R/2L,



Fig. 4R/4L
Columbia

Columbia

-64-134

Periodontia Instruments

Fig- 4R/4L,



Fig. 13/14
Columbia

Columbia

-64-135

Periodontia Instruments

Fig- 13/14,



Fig. 1/2
Langer

Langer

-64-136

Periodontia Instruments

Fig- 1/2,

Periodontia Instruments



Fig. 3/4
Langer

Langer

-64-137

Periodontia Instruments

Fig- 3/4,



Fig. 1/2
Barnhart

Barnhart

-64-138

Periodontia Instruments

Fig- 1/2,



Fig. 5/6
Barnhart

Barnhart

-64-139

Periodontia Instruments

Fig- 5/6,



Fig. ME-1
(5/11)

XU-64-140

Periodontia Instruments

Fig- ME-1,



Fig. 3/4
Hure (23/24)

XU-64-141

Periodontia Instruments

Fig- 3/4,



Fig. 1/2
Langer

Langer

-64-142

Periodontia Instruments

Fig- 1/2,

Periodontia Instruments



Fig. 3/4
Langer

Langer

-64-143

Periodontia Instruments

Fig- 3/4,



Fig. 1/2
Barnhart

Barnhart

-64-144

Periodontia Instruments

Fig- 1/2,



Fig. 5/8
Barnhart

Barnhart

-64-145

Periodontia Instruments

Fig- 5/6,



Fig. ME-1
[10/11]

XU-64-146

Periodontia Instruments

Fig- ME-1,



Fig. 3/4
Rive [23/24]

XU-64-147

Periodontia Instruments

Fig- 3/4,



Fig. 1/2
Gracey

Gracey

-64-148

Periodontia Instruments

Fig- 1/2,

Periodontia Instruments



Fig. 3/4
Gracey

Gracey

-64-149

Periodontia Instruments

Fig- 3/4,



Fig. 5/6
Gracey

Gracey

-64-150

Periodontia Instruments

Fig- 5/6,



Fig. 7/8
Gracey

Gracey

-64-151

Periodontia Instruments

Fig- 7/8,



Fig. 9/10
Gracey

Gracey

-64-152

Periodontia Instruments

Fig- 9/10,



Fig. 11/12
Gracey

Gracey

-64-153

Periodontia Instruments

Fig- 11/12,



Fig. 13/14
Gracey

Gracey

-64-154

Periodontia Instruments

Fig- 13/14,

Periodontia Instruments



Fig. 1/2
Gracey

Gracey

-64-155

Periodontia Instruments

Fig- 1/2,



Fig. 3/4
Gracey

Gracey

-64-156

Periodontia Instruments

Fig- 3/4,



Fig. 5/6
Gracey

Gracey

-64-157

Periodontia Instruments

Fig- 5/6,



Fig. 7/8
Gracey

Gracey

-64-158

Periodontia Instruments

Fig- 7/8,



Fig. 9/10
Gracey

Gracey

-64-159

Periodontia Instruments

Fig- 9/10,



Fig. 11/12
Gracey

Gracey

-64-160

Periodontia Instruments

Fig- 11/12,

Periodontia Instruments



Fig. 13/14
Gracey

Gracey

-64-161

Periodontia Instruments
Fig- 13/14,



Fig. 11/14
Gracey

Gracey

-64-162

Periodontia Instruments
Fig- 11/14,



Fig. 12/13
Gracey

Gracey

-64-163

Periodontia Instruments
Fig- 12/13,



Fig. 1
Ce Marco

Ce Marco

-64-164

Periodontia Instruments
Fig- 1,



Fig. 2
De Marco

De Marco

-64-165

Periodontia Instruments
Fig- 2,



Fig. 13/14
Anterior

Anterior

-64-166

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 17/18
Posterior

Posterior

-64-167

Periodontia Instruments
Fig- 17/18,



Fig. 11/14
Gracey

Gracey

-64-168

Periodontia Instruments
Fig- 11/14,



Fig. 12/13
Gracey

Gracey

-64-169

Periodontia Instruments
Fig- 12/13,



Fig. 1
De Marco

De Marco

-64-170

Periodontia Instruments
Fig- 1,



Fig. 2
De Marco

De Marco

-64-171

Periodontia Instruments
Fig- 2,



Fig. 13/14
Anterior

Anterior

-64-172

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 17/18
Posierior

**Posierior
-64-173**

Periodontia Instruments
Fig- 17/18,



Fig. 1
Goldman-Fox

**Goldman-Fox
XU-64-174**

Periodontia Instruments
Fig- 1,



Fig. 2
Goldman-Fox

**Goldman-Fox
XU-64-175**

Periodontia Instruments
Fig- 2,



Fig. 3
Goldman-Fox

**Goldman-Fox
XU-64-176**

Periodontia Instruments
Fig- 3,



Fig. 4
Goldman-Fox

**Goldman-Fox
XU-64-177**

Periodontia Instruments
Fig- 4,



Fig. 5
Goldman-Fox

**Goldman-Fox
XU-64-178**

Periodontia Instruments
Fig- 5,

Periodontia Instruments



Fig. 6
Goldman-Fox

Goldman-Fox
XU-64-179
Periodontia Instruments
Fig- 6,



Fig. 21
Goldman-Fox

Goldman-Fox
XU-64-180
Periodontia Instruments
Fig- 21,



Fig. 1
Goldman-Fox

Goldman-Fox
XU-64-181
Periodontia Instruments
Fig- 1,



Fig. 2
Goldman-Fox

Goldman-Fox
XU-64-182
Periodontia Instruments
Fig- 2,



Fig. 3
Goldman-Fox

Goldman-Fox
XU-64-183
Periodontia Instruments
Fig- 3,



Fig. 4
Goldman-Fox

Goldman-Fox
XU-64-184
Periodontia Instruments
Fig- 4,

Periodontia Instruments



Fig. 5
Goldman-Fox

Goldman-Fox
XU-64-185
Periodontia Instruments
Fig- 5,



Fig. 6
Goldman-Fox

Goldman-Fox
XU-64-186
Periodontia Instruments
Fig- 6,



Fig. 21
Goldman-Fox

Goldman-Fox
XU-64-187
Periodontia Instruments
Fig- 21,



Fig. 11/12
Mc Call

Mc Call
-64-188
Periodontia Instruments
Fig- 11/12,



Fig. 11A/12A
Mc Call

Mc Call
-64-189
Periodontia Instruments
Fig- 11A/12A,



Fig. 13/14
Mc Call

Mc Call
-64-190
Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 13S/14S
Mc Call

Mc Call

-64-191

Periodontia Instruments
Fig- 13S/14S,



Fig. 17S/18S
Mc Call

Mc Call

-64-192

Periodontia Instruments
Fig- 17S/18S,



Fig. 19/20
Mc Call

Mc Call

-64-193

Periodontia Instruments
Fig- 19/20,



Fig. 11/12
Mc Call

Mc Call

-64-194

Periodontia Instruments
Fig- 11/12,



Fig. 11A/12A
Mc Call

Mc Call

-64-195

Periodontia Instruments
Fig- 11A/12A,



Fig. 13/14
Mc Call

Mc Call

-64-196

Periodontia Instruments
Fig- 13/14,

Periodontia Instruments



Fig. 13S/14S
Mc Call

Mc Call

-64-197

Periodontia Instruments
Fig- 13S/14S,



Fig. 17S/18S
Mc Call

Mc Call

-64-198

Periodontia Instruments
Fig- 17S/18S,



Fig. 19/20
Mc Call

Mc Call

-64-199

Periodontia Instruments
Fig- 19/20,



Fig. 2/3
Taylor

Taylor

-64-200

Periodontia Instruments
Fig- 2/3,



Fig. 204

XU-64-201

Periodontia Instruments
Fig- 204,



Fig. 204S

XU-64-202

Periodontia Instruments
Fig- 204s,

Periodontia Instruments



Fig. 204S

XU-64-203
Periodontia Instruments
Fig- 204S,



Fig. 6
Crane-Kaplan

Crane-Kaplan
XU-64-204
Periodontia Instruments
Fig- 6,



Fig. 7/8
Younger-Good

Younger-Good
XU-64-205
Periodontia Instruments
Fig- 7/8,



Fig. 2/3
Taylor

Taylor
-64-206
Periodontia Instruments
Fig- 2/3,



Fig. 204

XU-64-207
Periodontia Instruments
Fig- 204,



Fig. 204S

XU-64-208
Periodontia Instruments
Fig- 204S,

Periodontia Instruments



Fig 204S

XU-64-209
Periodontia Instruments
Fig- 204S,



Fig. 6
Crane-Kaplan

Crane-Kaplan
XU-64-210
Periodontia Instruments
Fig- 6,



Fig 7/8
Younger-Good

Younger-Good
XU-64-211
Periodontia Instruments
Fig- 7/8,



Fig. 3/7
Hirschfeld

Hirschfeld
XU-64-212
Periodontia Instruments
Fig- 3/7,



Fig 30/33
Jaquette

Jaquette
-64-213
Periodontia Instruments
Fig- 30/33,



Fig. 31/32
Jaquette

Jaquette
-64-214
Periodontia Instruments
Fig- 31/32,

Periodontia Instruments



Fig. 31/32
Jaquette

Jaquette

-64-215

Periodontia Instruments
Fig- 31/32,



Fig. 3/7
Hirschfeld

Hirschfeld

XU-64-216

Periodontia Instruments
Fig- 3/7,



Fig. 30/33
Jaquette

Jaquette

-64-217

Periodontia Instruments
Fig- 30/33,



Fig. 31/32
Jaquette

Jaquette

-64-218

Periodontia Instruments
Fig- 31/32,



Fig. 34/35
Jaquette

Jaquette

-64-219

Periodontia Instruments
Fig- 34/35,



Fig. 37/38
Jaquette

Jaquette

-64-220

Periodontia Instruments
Fig- 37/38,

Periodontia Instruments



Fig. 1
Ochsenbein

Ochsenbein
XU-64-221
Periodontia Instruments
Fig- 1,



Fig. 2
Ochsenbein

Ochsenbein
XU-64-222
Periodontia Instruments
Fig- 2,



Fig. 1
Gardner

Gardner
-64-223
Periodontia Instruments
Fig- 1,



Fig. 2
Gardner

Gardner
-64-224
Periodontia Instruments
Fig- 2,



Fig. 5
Gardner

Gardner
-64-225
Periodontia Instruments
Fig- 5,



Fig. 5R
Gardner

Gardner
-64-226
Periodontia Instruments
Fig- 5R,

Periodontia Instruments



Fig. 6
Gardner

Gardner

-64-227

Periodontia Instruments

Fig- 6,



Fig. 1W
Wakef Old

Wakef Old

XU-64-228

Periodontia Instruments

Fig- 1W,



Fig. 2W
Wakef Old

Wakef Old

XU-64-229

Periodontia Instruments

Fig- 2W,



Fig. 3
Buckley

Buckley

-64-230

Periodontia Instruments

Fig- 3,



Fig. M9
Molt

Molt

-64-231

Periodontia Instruments

Fig- M9,



Fig. PR3
Prichard

Prichard

-64-232

Periodontia Instruments

Fig- PR3,

Periodontia Instruments



Fig. 1
Hopkins

Hopkins

-64-233

Periodontia Instruments

Fig- 1,



Fig. 2
Hopkins

Hopkins

-64-234

Periodontia Instruments

Fig- 2,



Fig. 5
Frerr

Frerr

-64-235

Periodontia Instruments

Fig- 15,



Fig. 16
Frerr

Frerr

-64-236

Periodontia Instruments

Fig- 16,



Fig. 3
Mead

Mead

-64-237

Periodontia Instruments

Fig- 3,



Fig. 23
Seldin retractor

Seldin retractor

XU-64-238

Periodontia Instruments

Fig- 23,

Periodontia Instruments



Williger
cm. 17

Williger

-64-239

Periodontia Instruments
cm 17,



Williger
cm. 16

Williger

-64-240

Periodontia Instruments
cm 16,



Williger
cm. 16

Williger

-64-241

Periodontia Instruments
cm 16,



Williger
cm. 13

Williger

-64-242

Periodontia Instruments
cm 13,



Fig 1S/2S
Sugarman

Sugarman

XU-64-243

Periodontia Instruments
Fig- 1S/2S,



Fig 3S/4S
Sugarman

Sugarman

XU-64-244

Periodontia Instruments
Fig- 3S/4S,

Periodontia Instruments



Fig. 1/2
TG Chisel

TG Chisel

-64-245

Periodontia Instruments

Fig- 1/2,



Fig. 1/2
TG Chisel

TG Chisel

-64-246

Periodontia Instruments

Fig- 1/2,



Fig. 1
TG Chisel

TG Chisel

-64-247

Periodontia Instruments

Fig- 1,



Fig. 2
TG Chisel

TG Chisel

-64-248

Periodontia Instruments

Fig- 2,



Fig. 1
(196)
Kirkland

Kirkland

-64-249

Periodontia Instruments

Fig- 1,



Fig. 2
(194)
Kirkland

Kirkland

-64-250

Periodontia Instruments

Fig- 2,

Periodontia Instruments



Fig. 1
Orban

Orban

-64-251

Periodontia Instruments
Fig- 1,




Fig. 2
Orban

Orban

-64-252

Periodontia Instruments
Fig- 2,

Filling instruments




Westcott

-65-101

Filling instruments

Fig- 1,




Westcott

-65-110

Filling instruments


Fig- 1,



XU-65-119

Filling instruments

Fig- 1, Fig- 2, Fig- 3, Fig- 4, Fig- 5,



XU-65-125

Filling instruments

Fig- 1, Fig- 2, Fig- 3, Fig- 1L, Fig- 2L, Fig- 3L,

Cement Spatulas



XU-66-101
No. 22 Cement Spatula



XU-66-102
No. 24A Cement Spatula, Rigid



XU-66-103
No. 24 Cement Spatula, Flexible



XU-66-104
No. 313 Cement Spatula,




XU-66-105
No. 5 Cement Spatula




XU-66-106
No. 7 Cement Spatula

Single/Double End Explorers



XU-67-101

Single Ended/ Double Ended Explorers
Fig- UT 1/2,



XU-67-112

Single Ended/ Double Ended Explorers
Fig- UT 1/2,

Mirror Handles



SIMPLE STEM PLANE 4

XU-68-113
MOUTH MIRROR
22 mm,



SIMPLE STEM PLANE 5

XU-68-114
MOUTH MIRROR
24 mm,



SIMPLE STEM PLANE 4

XU-68-115
MOUTH MIRROR
22 mm,



SIMPLE STEM PLANE 4

XU-68-116
MOUTH MIRROR
22 mm,



SIMPLE STEM PLANE 4

XU-68-117
MOUTH MIRROR
22 mm,



SECO / CONE-SOCKET

XU-68-118
Mirror Handles
125 mm,5,

Mirror Handles



OCTOGONAL CONE SOCKET

XU-68-119

Mirror Handles
120 mm, 4,



MARTIN ERGO

XU-68-120

Mirror Handles
135 mm, 5?,



Mirror Handles

XU-68-121

Mirror Handles
135 mm, 5?,



Mirror Handles

XU-68-122

Mirror Handle
135 mm, 5?,



Mirror Handles

XU-68-123

Mirror Handles
135 mm, 5?,



Mirror Handles

XU-68-124

Mirror Handles
135 mm, 5?,

Mirror Handles



Laryngeal Mirrors
XU-68-128
MIRROR HANDLES
7 mm 30°, Ø 4,



Laryngeal Mirrors
XU-68-129
MIRROR HANDLES
7 mm 30°, Ø 4,



Cotton Wool Container
XU-68-130
Cotton Wool Container
Ø 50 mm,



Napkin Holder
XU-68-131
Napkin Holder
Ø 50 mm,



MIRROR HANDLE CONE SOCKET
XU-53-721-01
MIRROR HANDLES
14 cm /5,



MIRROR HANDLE CONE SOCKET
XU-53-721-02
MIRROR HANDLES
14 cm /5,

Mirror Handles



MIRROR HANDLE CONE SOCKET

XU-53-723-01

MIRROR HANDLES
13 cm /5,



MIRROR HANDLE

XU-53-728-02

MIRROR HANDLES
12 cm /4,



MARTIN-GRIP

XU-53-746-02

12 cm /4?,



MIRROR HANDLE

XU-53-755-02

MIRROR HANDLES
14 cm /5 Ø 11,



MIRROR HANDLE

XU-53-757-08

MIRROR HANDLES
14 cm /5,



MIRROR

XU-53-702-

22

Mirror Handles



MIRROR
XU-53-703-22



RHODIUM MIRROR
XU-53-704-10
MIRROR



RHODIUM MIRROR
XU-53-704-18
MIRROR



CONE SOCKET
XU-53-706-22
MIRROR



CONE SOCKET
XU-53-707-22
MIRROR



MICRO-MOUTH MIRROR
XU-53-710-03
Fig. 1,


Mirror Handles



MICRO-MOUTH MIRROR

XU-53-710-05

Fig. 2,



RHODIUM MICRO MIRROR

XU-53-711-03

Fig.1 ? 3,0 mm,



HILGER

XU-53-715-

00

18 cm /7,



HILGER

-53-715-01

CHEEK RETRACTORS

Size 1,

Endodontic Instruments



Fig. 1
Spreader

Root Spreader fig 1
XU-69-101
Endodontic Instruments
Fig. 1,



Fig. 2
Spreader

Root Spreader fig 2
XU-69-102
Endodontic Instruments
Fig- 2,



Fig. 3
Spreader

Root Spreader Fig 3
XU-69-103
Endodontic Instruments
Fig- 3,



Fig. 4
Spreader

Root Spreader Fig 4
XU-69-104
Endodontic Instruments
Fig- 4,



Fig. 5
Spreader

Root Spreader Fig 5
XU-69-105
Endodontic Instruments
Fig- 5,



Fig. 6
Spreader

Root Spreader Fig 6
XU-69-106
Endodontic Instruments
Fig- 6,

Endodontic Instruments



Fig. 7
Spreader

Root Spreader Fig 7
XU-69-107
Endodontic Instruments
Fig- 7,



Fig. 1
Luks

Luks Fig 1
XU-69-108
Endodontic Instruments
Fig- 1,



Fig. 2

LUKS Fig 2
XU-69-109
Endodontic Instruments
Fig- 2,



Fig. 3

LUKS Fig 3
XU-69-110
Endodontic Instruments
Fig- 3,



Fig. 4

LUKS Fig 4
XU-69-111
Endodontic Instruments
Fig- 4,



Fig. 1

Vertical Condenser Fig 1
XU-69-112
Endodontic Instruments
Fig- 1,

Endodontic Instruments



Fig. 2

Vertical Condenser Fig 2

XU-69-113

Endodontic Instruments
Fig- 2,



Fig. 3

Vertical Condenser Fig 3

XU-69-114

Endodontic Instruments
Fig- 3,



Fig. 6/1

Root Canal Explorer Fig 6/1

XU-69-116

Endodontic Instruments
Fig- 6/1,



Fig. 6/17

Root Canal Explorer Fig 6/17

XU-69-117

Endodontic Instruments
Fig- 6/17,



Fig. 6/23

Root Canal Explorer Fig 6/23

XU-69-118

Endodontic Instruments
Fig- 6/23,



Fig. 31L

Root Canal Excavators Fig 31L

XU-69-119

Endodontic Instruments
Fig- 31L,

Endodontic Instruments



Fig. 32L

Root Canal Excavators Fig 32L

XU-69-120

Endodontic Instruments
Fig- 32L,



Fig. 33L

Root Canal Excavators Fig 33L

XU-69-121

Endodontic Instruments
Fig- 33L,



Fig. 6/1

Root Canal Explorers Fig 6/1

XU-69-122

Endodontic Instruments
Fig- 6/1,



Fig. 6/17

Root Canal Explorers Fig 6/17

XU-69-123

Endodontic Instruments
Fig- 6/17,



Fig. 6/23

Root Canal Explorers Fig 6/23

XU-69-124

Endodontic Instruments
Fig- 6/23,



Fig. 31L

Root Canal Excavators Fig 31L

XU-69-125

Endodontic Instruments
Fig- 31L,

Endodontic Instruments



Fig. 32L

Root Canal Excavators Fig 32L

XU-69-126

Endodontic Instruments
Fig- 32L,



Fig. 33L

Root Canal Excavators Fig 33L

XU-69-127

Endodontic Instruments
Fig- 33L,

Separating Spatulas



Heidemann
XU-70-101
Separating Spatulas
Fig- 1,



Heidemann
XU-70-102
Separating Spatulas
Fig- 2,



Heidemann
XU-70-104
Separating Spatulas
Fig- 4,



Heidemann
XU-70-105
Separating Spatulas
Fig- 5,



Heidemann
XU-70-106
Separating Spatulas
Fig- 6,



Heidemann
XU-70-107
Separating Spatulas
Fig- 7,

Separating Spatulas



Heidemann
XU-70-108
Separating Spatulas
Fig- 8,



Heidemann
XU-70-109
Separating Spatulas
Fig- 9,



Heidemann
XU-70-110
Separating Spatulas
Fig- 1,



Heidemann
XU-70-111
Separating Spatulas
Fig- 2,



Heidemann
XU-70-112
Separating Spatulas
Fig- 3,

Excavators



Fig. 0

Excavators Fig 0

XU-71-101
Excavators
Fig- 0,



Fig. 1

Excavators Fig 1

XU-71-102
Excavators
Fig- 1,



Fig. 2

Excavators Fig 2

XU-71-103
Excavators
Fig- 2,



Fig. 3

Excavators Fig 3

XU-71-104
Excavators
Fig- 3,



Fig. 4

Excavators Fig 4

XU-71-105
Excavators
Fig- 4,



Fig. 5

Excavators Fig 5

XU-71-106
Excavators
Fig- 5,

Excavators



Fig. 6

Excavators Fig 6

XU-71-107
Excavators
Fig- 6,



Fig. 0

Excavators Fig 0

XU-71-108
Excavators
Fig- 0,



Fig. 1

Excavators Fig 1

XU-71-109
Excavators
Fig- 1,



Fig. 2

Excavators Fig 2

XU-71-110
Excavators
Fig- 2,



Fig. 3

Excavators Fig 3

XU-71-111
Excavators
Fig- 3,



Fig. 4

Excavators Fig 4

XU-71-112
Excavators
Fig- 4,

Excavators



Fig. 5

Excavators Fig 5

XU-71-113
Excavators
Fig- 5,



Fig. 6

Excavators Fig 6

XU-71-114
Excavators
Fig- 6,



Fig. 6

Excavators Fig 6

XU-71-115
Excavators
Fig- 6,



Fig. 14

Excavators Fig 14

XU-71-116
Excavators
Fig- 14,



Fig. 17

Excavators Fig 17

XU-71-117
Excavators
Fig- 17,



Fig. 17W

Excavators Fig 17W

XU-71-118
Excavators
Fig- 17W,

Excavators



Fig. 18

Excavators Fig 18

XU-71-119
Excavators
Fig- 18,



Fig. 18W
Oval spoon

Oval spoon Fig 18W

XU-71-120
Excavators
Fig- 18W,



Fig. 19

Excavators Fig 19

XU-71-121
Excavators
Fig- 19,



Fig. 19W
Oval spoon

Oval spoon Fig 19W

XU-71-122
Excavators
Fig- 19W,



Fig. 6

Excavators Fig 6

XU-71-123
Excavators
Fig- 6,



Fig. 14

Excavators Fig 14

XU-71-124
Excavators
Fig- 14,

Excavators



Fig. 17

Excavators Fig 17

XU-71-125
Excavators
Fig- 17,



Fig. 17W
Oval spoon

Oval spoon Fig 17W

XU-71-126
Excavators
Fig- 17W,



Fig. 18W
Oval spoon

Oval spoon 18W

XU-71-127
Excavators
Fig- 18W,



Fig. 19

Excavators Fig 19

XU-71-129
Excavators
Fig- 19,



Fig. 19W
Oval spoon

Oval spoon Fig 19W

XU-71-130
Excavators
Fig- 19W,



Fig. 38/39

Excavators Fig 38/39

XU-71-131
Excavators
Fig- 38/39,

Excavators



Fig. 63/64

Excavators Fig 63/64

XU-71-133
Excavators
Fig- 63/64,



Fig. 65/66

Excavators Fig 65/66

XU-71-134
Excavators
Fig- 65/66,



Fig. 31L

Excavators Fig 31L

XU-71-135
Excavators
Fig- 31L,



Fig. 32L

Excavators Fig 32L

XU-71-136
Excavators
Fig- 32L,



Cleoid / Discoid

XU-71-137
No. 3/6 Excavator, double end



Cleoid / Discoid

XU-71-138
No. 89/92 Carver - Discoid \"/> Cleoid,
octagonal handle

Scalers



Fig. U15
Towner

Towner Fig U15

XU-72-104

Scalers
Fig- U15,



Fig. U15/33
Towner Jaquette

Towner Jaquette Fig- U15/33

XU-72-105

Scalers
Fig- U15/33,



Fig. U15/30
Towner Jaquette

Towner Jaquette Fig- 15/30

XU-72-106

Scalers
Fig- 15/30,



Fig. H5/33

Fig- H5/33

XU-72-107

Scalers
Fig- H5/33,



Fig. H6/H7

Fig- H6/H7

XU-72-108

Scalers
Fig- H6/H7,



Fig. 1/2

Fig- 1/2

XU-72-

109

Scalers
Fig- 1/2,

Scalers



Fig. U15/33
Towner Jaquette

Towner Jaquette Fig- U15/33

XU-72-110

Scalers
Fig- U15/33,



Fig. J15/30
Towner Jaquette

Towner Jaquette Fig- U15/30

XU-72-111

Scalers
Fig- U15/30,



Fig. 1

Fig- 1

XU-72-

112

Scalers
Fig- 1,



Fig. 33

XU-72-117

Scalers

Fig- 33, Fig- 31, Fig- 32, Fig- 34, Fig- 35,



Fig. H5
Hygienist

Hygienist Fig- H5

XU-72-122

Scalers

Fig- H5,

Wax & Modelling Carvers



P.K. THOMAS CARVERS FIG. 1

XU-73-101

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.1



P.K. THOMAS CARVERS FIG. 2

XU-73-102

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.2



P.K. THOMAS CARVERS FIG. 3

XU-73-103

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.3



P.K. THOMAS CARVERS FIG. 4

XU-73-104

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.4



P.K. THOMAS CARVERS FIG. 5

XU-73-105

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.5



P.K.THOMAS INSTRUMENTS

XU-73-106

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.5 SILVER

Wax & Modelling Carvers



P.K.THOMAS INSTRUMENTS

XU-73-107

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.4 RED



P.K.THOMAS INSTRUMENTS

XU-73-108

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.3 BLUE



P.K.THOMAS INSTRUMENTS

XU-73-109

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.2 GOLD



P.K.THOMAS INSTRUMENTS

XU-73-110

P.K.THOMAS WAX AND MODELLING
CARVERS FIG.1 GREEN



XU-73-111

FAHNENSTOCK Wax Carver \\\\\\\\\\\\\\\\" Knife,
5\\\\\\\\\\\\\\\\\\\", small



XU-73-112

FAHNENSTOCK Wax Carver \\\\\\\\\\\\\\\\" Knife, 6-
3/4\\\\\\\\\\\\\\\\\\\", large

Impression Trays



Impression Trays

XU-74-101

Impression Trays
66x52 mm,



Impression Trays

XU-74-102

Impression Trays
62x50 mm,



Impression Trays

XU-74-103

Impression Trays
64x50 mm,



Impression Trays

XU-74-104

Impression Trays
65x46 mm,



Impression Trays

XU-74-105

Impression Trays
69x54 mm,



Impression Trays

XU-74-106

Impression Trays
65x46 mm,

Impression Trays



Impression Trays

XU-74-107

Impression Trays
70x48 mm,



Impression Trays

XU-74-108

Impression Trays
70x55 mm,



Impression Trays

XU-74-109

Impression Trays
59x42 mm,



Impression Trays

XU-74-110

Impression Trays
55x37 mm,



Impression Trays

XU-74-111

Impression Trays
61x42 mm,



Impression Trays

XU-74-112

Impression Trays
44 mm,

Impression Trays



Impression Trays

XU-74-113

Impression Trays
44 mm,



Impression Trays

XU-74-114

Impression Trays
44 mm,



Impression Trays

XU-74-115

Impression Trays
66x52 mm,



Impression Trays

XU-74-116

Impression Trays
65x46 mm,



Impression Trays

XU-74-117

Impression Trays
59x42 mm,



Impression Trays

XU-74-118

Impression Trays
69x54 mm,

Impression Trays



Impression Trays

XU-74-119

Impression Trays
70x48 mm,



Impression Trays

XU-74-120

Impression Trays
70x55 mm,



Impression Trays

XU-74-121

Impression Trays
55x37 mm,



Impression Trays

XU-74-122

Impression Trays
61x42 mm,



Impression Trays

XU-74-123

Impression Trays
61x42 mm,



Impression Trays

XU-74-124

Impression Trays
61x42 mm,

Impression Trays



Impression Trays

XU-74-125

Impression Trays

Dental Syringes



XU-75-101
Dental Syringes
55 x 37 mm,



XU-75-102
Dental Syringes
55 x 37 mm,




Self Aspirating 1.8cc
XU-75-103
Dental Syringes
55 x 37 mm,



Cylinder Cartridge Syringe
XU-75-104
Dental Syringes
55 x 37 mm,



XU-75-105
Dental Syringes
55 x 37 mm,



XU-75-106
Dental Syringes
55 x 37 mm,

Dental Syringes



XU-75-107
Dental Syringes
55 x 37 mm,



C-W 1.8cc
XU-75-108
Dental Syringes
55 x 37 mm,



XU-75-109
Dental Syringes
55 x 37 mm,



XU-75-110
Dental Syringes
55 x 37 mm,



XU-75-111
Dental Syringes
55 x 37 mm,



XU-75-112
Dental Syringes
55 x 37 mm,

Dental Syringes



XU-75-113
Dental Syringes
55 x 37 mm,



XU-75-118
Dental Syringes
Fig .1, Fig .2, Fig .3,



Dental Syringes
XU-75-119
Aspira-Plus
1.8 ml,



Dental Syringes
XU-75-120
Dental Syringes
1.8 ml,



Dental Syringes
XU-75-122
Dental Syringes
1,8 ml,



Dental Syringes
XU-75-124
Dental Syringes
1,8 ml,

Dental Syringes



Dental Syringes
XU-75-126
Dental Syringes
1,8 ml,



Dental Syringes
XU-75-127
Dental Syringes
1,8 ml,



Dental Syringe
XU-75-128
Aspira-Plus
1.8 ml,



Dental Syringe
XU-75-129



Dental Syringe
XU-75-130




Dental Syringe
XU-75-131

Dental Syringes



Dental Syringe
XU-75-132
1 ml,

Cavity Preparation




Wedelstaet

XU-76-101

Cavity Preparation

Fig- 1, Fig- 2, Fig- 3, Fig- 4, Fig- 1/2, Fig- 3/4, Fig- 5/6,



Wedelstaet

XU-76-112

Cavity Preparation

Fig- 1, Fig- 2, Fig- 3, Fig- 4, Fig- 1/2, Fig- 3/4, Fig- 5/6,

Explorers & Probes



11 Briault D/E Explorer
XU-77-101

FIG. # 11 BRIAULT D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



11/12 Odu D/E Explorer
XU-77-102

FIG. # 11/12 OLD DOMINION UNIVERSIT D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



11/12 Odu D/E After Five Explorer
XU-77-103

FIG. # 11/12 ODU D/E AFTER FIVE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2A Pigtail D/E Explorer
XU-77-104

FIG. # 2A PIGTAIL D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2H D/E Explorer
XU-77-105

FIG. # 2H DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



2R/2L D/E Explorer
XU-77-106

FIG. # 2R/2L D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE

Explorers & Probes



3 Pigtail D/E Explorer
XU-77-107

FIG. # 3 PIGTAIL D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3CH Cowhorn D/E Explorer
XU-77-108

FIG. # 3CH COWHORN D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



3ES D/E Explorer
XU-77-109

FIG. # 3ES DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



5 D/E Explorer
XU-77-110

FIG. # 5 D/E EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



6/3 D/E Explorer
XU-77-111

FIG. # 6/3 DE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE



17 S/E Explorer
XU-77-112

FIG. # 17 SE EXPLORER. MADE OF SOLID STAINLESS STEEL HANDLE

Explorers & Probes



23 Shepherd's Hook Explorer
XU-77-113

FIG. # 23 SHEPHERD'S HOOK S/E
EXPLORER. MADE OF SOLID STAINLESS
STEEL HANDLE



3 S/E Explorer
XU-77-114

FIG. # 3 SE EXPLORER. MADE OF SOLID
STAINLESS STEEL HANDLE



3A S/E Explorer
XU-77-115

FIG. # 3A S/E EXPLORER. MADE OF
SOLID STAINLESS STEEL HANDLE



54 S/E Explorer
XU-77-116

FIG. # 54 S/E EXPLORER. MADE OF
SOLID STAINLESS STEEL HANDLE

Crown Removers



XU-78-101

A tips for MORREL Crown Remover



XU-78-102

B tips for MORREL Crown Remover



XU-78-103

C tips for MORREL Crown Remover



XU-78-104

MORREL Crown Remover, with 3 points



XU-78-105

WYNMAN Carbide Crown Gripper, 5-1/2",
without ratchet, with TC jaws

X-Rays Holders



X-Ray Film Hangers 6 Clips

XU-79-101

X-RAY FILM HANGERS 6 CLIPS MADE OF
STAINLESS STEEL



X-Ray Film Hangers 10 Clips

XU-79-102

X-RAY FILM HANGERS 10 CLIPS MADE
OF STAINLESS STEEL

Matrix Retainers



Matrix Retainer Toffilemire

XU-80-101

MATRIX BAND RETAINER TOFFILEMIRE.
AVAILABLE IN STRAIGHT & CONTRA.
STAINLESS STEEL



**Toffilemire Matrix Retainer XUnior /
Child**

XU-80-102

MATRIX BAND RETAINER TOFFILEMI
XUNIOR / CHILD. AVAILABLE IN
STRAIGHT & CONTRA. STAINLESS
STEEL



Toffilemire Matrix Retainer Senior / Adult

XU-80-103

MATRIX BAND RETAINER TOFFILEMIRE
SENIOR / CHILD. AVAILABLE IN
STRAIGHT & CONTRA. STAINLESS
STEEL



Toffilemire Matrix Retainer Universal

XU-80-104

MATRIX BAND RETAINER TOFFILEMIR
UNIVERSAL. AVAILABLE IN STRAIGHT &
CONTRA. STAINLESS STEEL



Matrix Band Retainer Ivory 8

XU-80-105

MATRIX BAND RETAINER IVORY 8,
AVAILABLE IN 5MM, 6MM, 7MM.
STAINLESS STEEL



Matrix Band Retainer Ivory Plain

XU-80-106

MATRIX RETAINER IVORY PLAIN /
WITHOUT LATERAL SHIELD (JAW).
STAINLESS STEEL

Matrix Retainers



Matrix Retainer Ivory With Shield

XU-80-107

MATRIX RETAINER IVORY WITH LATERAL
SHIELD (JAW). STAINLESS STEEL



Matrix Retainer Cheek Protector

XU-80-108

MATRIX BAND RETAINER CHEECK
PROTECTOR. STAINLESS STEEL



Matrix Band Retainer Nystrom Left

XU-80-109

MATRIX BAND RETAINER NYSTROM
LEFT. STAINLESS STEEL



Matrix Band Retainer Nystrom Right

XU-80-110

MATRIX BAND RETAINER NYSTROM
RIGHT. STAINLESS STEEL



**Matrix Band Retainer Siqveland (Ash
Narrow)**

XU-80-111

MATRIX BAND RETAINER SIQVELAND
(ASH NARROW). STAINLESS STEEL



**Matrix Band Retainer Siqveland (Ash
Wide)**

XU-80-112

MATRIX BAND RETAINER SIQVELAND
(ASH WIDE). STAINLESS STEEL

Matrix Retainers



**Matrix Band Retainer Ivory Didner
(AdXUstable)**

XU-80-113

MATRIX BAND RETAINER IVORY DIDNER
(AdXUstable)



Matrix Retainer Ivory Separator

XU-80-114

MATRIX BAND RETAINER IVORY
SEPARATOR. STAINLESS STEEL



Matrix Retainer Elliott Separator

XU-80-115

MATRIX RETAINER ELLIOTE
SEPARATOR. STAINLESS STEEL



Matrix Retainer Ivory 8 N

XU-80-116

MATRIX BAND RETAINER IVORY 8N.
MADE OF STAINLESS STEEL




Matrix Retainer Didner 8 P

XU-80-117

MATRIX BAND RETAINER IVORY 8P.
MADE OF STAINLESS STEEL

Articulators



**PLASTERLESS ARTICULATOR GALETTI
(WHITE)
XU-81-101**

HIGH TECH PLASTERLESS ARTICULATOR GALETTI WITH ADJUSTABLE SCREWS A REARS & FRONT. SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING. ALL ANATOMICAL MOVEMENTS POSSIBLE



**Vibratory Polished Articulator Galetti
XU-81-102**

GALETTI ARTICULATOR IS SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING. EASY TO USE, ALLOWING QUICK FIXATION OF MODEL. IMMEDIATE ADJUSTMENT OF THE OCCLUSAL & CORRECTS IMITATION OF CENTRAL RELATION. ALL ANATOMICAL MOVEMENTS POSSIBLE. VIBRATORY POLISHED



**ANATOMICAL ARTICULATOR (ASA)
XU-81-103**

ANATOMICAL ARTICULATOR WITH BEARING ROD & INCISIVE PIN, 4 MOUNTING PLATES & SUPPORTING BACK-PIN FOR 45° VIEW, OPENING 125°, INTERCONDYLAR DISTANCE 110 MM - SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING



**UNIVERSAL ARTICULATOR
XU-81-104**

ANATOMICAL ARTICULATOR, MADE OF ALUMINUM PRESSOFUSION OF ONLY 590/GR OF WEIGHT. REALIZED FOR EASY VISIBILITY TO THE LINGUAL SURFACE & FOR A PERFECT STABILITY ALSO IN RESERVED POSITION. THE RETURN TO THE HINGES WITH RECALL SPRING IS AUTOMATIC. OPENING 125°, WITH 4 MOUNTING PLATES (AVAILABLE



**Gnatus Pro JP30 Articulator
XU-81-105**

The Gnatus Pro JP30 Articulator offers the best performance with an excellent degree of accuracy. Its functionality makes easier both the procedures of taking and transferring the patient's records and the process of mounting the models. Records can be taken and transferred accurately and fast. Light and resistant, makes operation easier. Its symmetrical and centered opening and



**ANATOMICAL ARTICULATOR
XU-81-106**

ANATOMICAL ARTICULATOR WITH INCISIVE PIN, 4 MOUNTING PLATES AND SUPPORTING BACK-PIN FOR 45° VIEW, OPENING 125°, INTERCONDYLAR DISTANCE 110 MM. - AVAILABLE IN ASSORTED COLOR

Articulators




KEYSTONE BLACK ARTICULATOR WIRED
XU-81-107

KEYSTONE DELUXE UNIVERSAL TYPE-A ARTICULATOR. INTERCHANGEABLE, INCISAL PIN INCLUDED, LATERAL EXCURSIONS.



KEYSTONE DELUXE MAGNET ARTICULATOR
XU-81-108

KEYSTONE DELUXE MAGNETIC DENTURE ARTICULATOR WITH PLATES, INTERCHANGEABLE, INCISAL PIN INCLUDED. MAGNETIC SPLIT CASTING ENSURES EASY & EXACT SETTING, LATERAL EXCURSIONS



PROFESSIONAL OCCULUSOR ARTICULATOR
XU-81-109

PROFESSIONAL OCCULUSOR, VIBRATOR POLISHED ARTICULATOR - ALUMINUM, BRASS, STAINLESS STEEL




Magicart-1 Articulator
XU-81-110

LIGHTWEIGHT, HIGH-QUALITY MAGIC ART ARTICULATOR SECURELY HOLDS QUADRANT MODELS, CONSERVE MATERIALS & TIME! QUICKLY & EASILY MOUNT PRE-CAST DENTAL MODELS WITH THE HIGH QUALITY. WITH MAGICART, YOU CAN SECURELY MOUNT & REMOVE ANY SIZE QUADRANT MODE CASTS WITH AN EXCLUSIVE SIDE GRIP. SERRATED C-CLAMPS LOCK MODELS TIGHTLY WITH A QUICK ADJUSTMENT &



FULL DENTURE ARTICULATOR (ROUND)
XU-81-111
FULL DENTURE CHROME PLATED ARTICULATOR (3 SLOT ROUND)



MODEL 'M' HI TOP STEPHAN ARTICULATOR
XU-81-112
FULL DENTURE CHROMIUM PLATED ARTICULATOR. PROVIDES ADDITIONAL VERTICAL DIMENSION FOR ACCEPTANCE OF HIGHER MODELS. 49

Articulators



Chromium Plated Plain Round

XU-81-113

Nickle Chromium Plated Plain Round
Articulator



UNIVERSAL ARTICULATOR (YELLOW)

XU-81-114

UNIVERSAL ARTICULATOR OPENING 125°
WITH 4 MOUNTING PLATES - CONDYLAR
INCLINATION FIXED AT 30°,
INTERCONDYLAR DISTANCE AT 110MM,
BENNET ANGLE FIXED AT 7°, DISTANCE
BETWEEN MOUNTING PLATES 85 MM.
SOLIDLY CONSTRUCTED OF METAL
ALLOY CASTING. Available In Assorted
Colors



FOSTER DENTURE ARTICULATOR

XU-81-115

THIS ARTICULATOR IS EQUIPPED WITH
AN ADJUSTABLE ANTERIOR PIN STOP
& POSTERIOR SCREW STOP, AND
CONSTRUCTED OF A SLIGHTLY HEAVY
GAUGE BRASS GIVING THE
ARTICULATOR A MORE SOLID FEEL.



FORK CHROME PLATED ARTICULATOR

XU-81-116

NICKLE CHROMIUM PLATED
ARTICULATOR FORK PATTERN



HINGE ARTICULATOR FULL ARCH

XU-81-117

AN ECONOMICAL & SIMPLE CHROMIUM
PLATED HINGE ARTICULATOR FOR
MOUNTING FULL ARCH CASES. LOW



**QUADRANT / CROWNED BRIDGE
ARTICULATOR**

XU-81-118

BRASS HINGE CROWNED BRIDGE
ARTICULATOR FOR MOUNTING
QUADRANT CASES. CROWN & BRIDGE
ARTICULATOR IS POPULAR FOSTER # 57
DESIGN, ELIMINATES THE NEED TO
BREAK OFF MODELS

Articulators



UNIVERSAL ARTICULATOR (BLUE)

XU-81-119

UNIVERSAL ARTICULATOR OPENING 125° WITH 4 MOUNTING PLATES - CONDYLAR INCLINATION FIXED AT 30°, INTERCONDYLAR DISTANCE AT 110MM, BENNET ANGLE FIXED AT 7°, DISTANCE BETWEEN MOUNTING PLATES 85 MM. SOLIDLY CONSTRUCTED OF METAL ALLOY CASTING.



FOSTER / APEX # 2 ARTICULATOR

XU-81-120

FOSTER / APEX # 2 ARTICULATOR, THE WIDTH AT THE JOINT IS 1-3/4" WITH A VERTICAL PARALLEL HEIGHT OF 2-1/8". THE SMALL SIZE OF THIS ARTICULATOR ALLOWS THE VERSATILITY TO USE FOR FULL OR PARTIAL WORK. THE FINISH IS BRIGHT CHROME. HIGH BACK



MORTIMER-I OCULADER ARTICULATOR

XU-81-121

PLAINLINE BALL BEARING MORTIMER-OCULADER ARTICULATOR. ALL FILLED WITH HIGH CARBON STEEL BALLS FOR SMOOTH SETTING & ARTICULATION OF DENTURE



MORTIMER-II OCULADER ARTICULATOR

XU-81-122

PLAINLINE BALL BEARING MORTIMER-OCULADER ARTICULATOR. ALL FILLED WITH HIGH CARBON STEEL BALLS. SEMI-ADJUSTABLE ARTICULATOR ALLOWS THE ESTABLISHMENT OF MORE ACCURATE OR CUSTOMISED LATERAL & PROTRUSIVE AS WELL AS CENTRIC OCCLUSION



OCCLUSOR, WIRE ARTICULATOR

XU-81-123

OCCLUSOR, WIRE ARTICULATOR - STAINLESS STEEL & ALUMINUM



ASA Articulator 5002

XU-81-124

ASA ARTICULATOR WITH LATERAL MOMENT. MADE OF ALUMINUM DIE CASTING WITH 2 MOUNTING PLATES. ALSO AVAILABLE IN VIBRATORY POLISHED PROCESS



Model 'A' Crescent Stephan Articulator

XU-81-125

CRESCENT STEPHAN MODEL "A" ARTICULATOR IS USED FOR PARTIALS, CENTER SPRING WITH ADJUSTABLE TENSION, SMOOTH BRASS CASTINGS WITH NICKEL PLATING, CENTER BAR IN BOTH UPPER & LOWER BOW FOR MAXIMUM PLASTER RETENTION & STRENGTH. WIDTH AT THE JOINT OF 29, & VERTICAL SPAN WHEN BOWS ARE PARALLEL OF 2 1/8". REPRODUCE ALL ANATOMICAL MOVEMENTS. UPPER AND LOWER SECTIONS MAY BE SEPARATE OR ASSEMBLED QUICKLY. OPEN CONSTRUCTION PROVIDES FREE ACCESS TO LINGUAL SURFACES. ANATOMICAL IN EVERY RESPECT



Bite Fork / Occlusal Alignment Plane

XU-81-126

BITE FORK IS A RELIABLE INSTRUMENT FOR CHECKING THE CORRECT POSITION OF THE MASTICATORY PLANE, ACHIEVING PARALLELISM OF THE OCCLUSAL PLANE WITH THE CAMPER'S LINE & THE BIPUPILLARY LINE IS FUNDAMENTAL TO THE CONSTRUCTION OF DENTURES

Miscellaneous Laboratory Instruments



XU-82-101

Bandage and Utility Scissors, 7-1/2" (19. cm), needle destroyer, serrated blade, autoclavable



XU-82-102

Plaster and Compound Knife



XU-82-103

Soldering Tweezers, curved



XU-82-104

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 3, narrow, fine, 4-3/4" (12.1 cm)



XU-82-105

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 5, Super fine, 4-3/8" (11.2 cm)



XU-82-106

SWISS Jeweler style Forceps, non-magnetic stainless steel, style 7, fine, curved, 4-1/2" (11.4 cm)

Miscellaneous Laboratory Instruments



Calipers
XU-82-107
Calipers



Calipers
XU-82-108
Calipers



Steel Ruler
XU-82-109
Steel Ruler



Forceps for articulation paper
XU-82-110
Forceps for articulation paper